

Analytical Letters

Analytical

Letters

ISSN: 0003-2719 (Print) 1532-236X (Online) Journal homepage: https://www.tandfonline.com/loi/lanl20

# Investigating the performance characteristics of the barrier discharge ionization detector and comparison to the flame ionization detector for the gas chromatographic analysis of volatile and semivolatile organic compounds

# Maria Antoniadou, George A. Zachariadis & Erwin Rosenberg

**To cite this article:** Maria Antoniadou, George A. Zachariadis & Erwin Rosenberg (2019) Investigating the performance characteristics of the barrier discharge ionization detector and comparison to the flame ionization detector for the gas chromatographic analysis of volatile and semivolatile organic compounds, Analytical Letters, 52:17, 2822-2839, DOI: 10.1080/00032719.2019.1628247

To link to this article: https://doi.org/10.1080/00032719.2019.1628247

9	© 2019 The Author(s). Published with license by Taylor and Francis Group, LLC	Published online: 03 Jul 2019.
	Submit your article to this journal $ arsigma^{\!$	Article views: 1571
Q	View related articles 🗷	View Crossmark data 🗹
ආ	Citing articles: 1 View citing articles 🖸	

#### GAS CHROMATOGRAPHY

OPEN ACCESS Check for updates

Taylor & Francis

Taylor & Francis Group

# Investigating the performance characteristics of the barrier discharge ionization detector and comparison to the flame ionization detector for the gas chromatographic analysis of volatile and semivolatile organic compounds

Maria Antoniadou<sup>a,b</sup> (b), George A. Zachariadis<sup>a</sup> (b), and Erwin Rosenberg<sup>b</sup> (b)

<sup>a</sup>Laboratory of Analytical Chemistry, Department of Chemistry, Aristotle University, Thessaloniki, Greece; <sup>b</sup>Institute of Chemical Technology and Analytics, Vienna University of Technology, Vienna, Austria

#### ABSTRACT

The response behavior and performance characteristics of the recently introduced barrier discharge ionization detector (BID) for gas chromatography (GC-BID) were investigated by analyzing different classes of organic compounds such as alcohols, alkanes, cycloaliphatic compounds, polycyclic aromatic hydrocarbons (PAHs), and others. The results obtained by GC-BID were compared with those of gas chromatography with flame ionization detection (GC-FID), aiming to demonstrate the particular merits of the new BID detector over the wellestablished FID. The response of the BID not only was found to be strongly dependent on the detector settings, but also shows a high dependence on the analyte class and the individual analyte. The sensitivity of the BID detector compared to the FID was higher by a factor of ca. 4 on average when considering all compounds analyzed. The relative standard deviation (RSD) was better than 5% for the majority of the cases. The BID detector showed better precision (lower RSD) in comparison with the FID for the investigated compounds. Linear calibrations were obtained for the analytes over more than four orders of magnitude with coefficients of determination typically higher than 0.999 and the limits of detection varied from 0.04 to 1.48 ng/s for the GC-BID.

#### **ARTICLE HISTORY**

Received 26 February 2019 Accepted 3 June 2019

#### **KEYWORDS**

Barrier discharge ionization detector (BID); flame ionization detector (FID); gas chromatography (GC); dielectric barrier discharge (DBD)

### Introduction

More than 75 years after its invention (Martin and Synge 1941), gas chromatography still maintains its role as the most important and most powerful analytical separation technique for volatile and semivolatile compounds. Its high separation efficiency, flexibility, and sensitivity are among the most commonly named assets of this technique. Much of the success of modern gas chromatography is due to the availability of a large number of sensitive detectors which, depending on the particular use, may be highly selective (such

© 2019 The Author(s). Published with license by Taylor and Francis Group, LLC.

CONTACT Erwin Rosenberg 🐼 erwin.rosenberg@tuwien.ac.at 💽 Institute of Chemical Technologies and Analytics, Vienna University of Technology, Getreidemarkt 9/164, A-1060 Vienna, Austria

Color versions of one or more of the figures in the article can be found online at www.tandfonline.com/lanl.

This paper is part of a special issue of papers presented at the Eleventh Aegean Analytical Chemistry Days (AACD2018), 25–29 September 2018, in Chania, Crete, Greece, organized by Prof. Antony C. Calokerinos, Prof. Nikos Thomaidis, and Prof. Nikos Lydakis-Simantiris.

This is an Open Access article distributed under the terms of the Creative Commons Attribution-NonCommercial-NoDerivatives License (http://creativecommons.org/licenses/by-nc-nd/4.0/), which permits non-commercial re-use, distribution, and reproduction in any medium, provided the original work is properly cited, and is not altered, transformed, or built upon in any way.

as the spectrometric detectors), or on the contrary universally applicable (such as the flame ionization detector, FID, or the thermal conductivity detector, TCD) (Rosenberg 2006).

One of the most recently introduced GC detectors is the so-called barrier-discharge ionization detector (BID). Being related to the FID through its detection principle, it also responds—in contrast to the former—to  $N_2$  and  $O_2$  and oxidized organic substances. The BID is based on a dielectric-barrier discharge (DBD) as the ionization source. The DBD is a low-temperature helium plasma, operated in a state of thermal nonequilibrium. It is produced by an electric discharge between two electrodes, at least one of which has a dielectric material on its surface (Liu et al. 2011). The discharge is initiated by a high voltage alternating current (AC) that is applied to the dielectric material. If the applied voltage surpasses the breakdown voltage of the plasma gas (also called the discharge gas, in this case He), a discharge is initiated between the two electrodes. The produced metastable species and photons ( $\Delta E$  up to 24.6 eV) are responsible for the ionization of the analytes. Because the electrode/dielectric barrier/electrode arrangement functions as a capacitor, each discharge is self-limiting and the formation of an electric arc between the electrodes (often associated with the sputtering of electrode material) is prevented (Gras et al. 2006; Jiang et al. 2014). The frequency and the voltage range from a few Hz to MHz and from 1 to 100 kV, respectively (Li et al. 2011a).

Originally described by Siemens (1857), this type of plasma was initially only adopted in production processes due to its low plasma temperature and power consumption, its simple manufacture and long operational use (Kogelschatz 2003; Hu et al. 2011). Only more recently, analytical and spectroscopic applications have become more abundant, and the DBD has been used in a number of analytical and spectroscopic applications (Guo et al. 2015; Schütz et al. 2015)

Typical applications of the DBD in spectroscopy are as an atomization (Straka et al. 2018; Kratzer et al. 2018), excitation (Brandt et al. 2016), and ionization source (Hagenhoff and Hayen 2018; Huba, Mirabelli, and Zenobi 2018). The DBD has also been incorporated as ionization source in some GC detectors, e.g., as GC photoionization detector on a chip (Zhu et al. 2016) and as an atmospheric pressure ionization detector providing detection in the low parts per billion range (Kirk, Last, and Zimmermann 2017). The combination of the DBD with miniaturized photodiode array detection has been reported particularly for element-specific detection in a number of applications (Li et al. 2011; Han et al. 2014; Li, Jiang, and Hou 2015; Jiang et al. 2015).

Also, the dielectric-barrier ionization detector that has recently become commercially available for GC-BID (Shinada et al. 2012) is based on the principle of dielectric barrier discharge ionization with He as a plasma gas. The detector, consisting of a plasma generator and the charge collector, are operated at medium frequency and high voltage (5–30 kHz, 5 to 10 kV) applied to a central ring electrode while two outer electrodes are grounded. The quartz tube acts as a dielectric barrier across which a discharge is initiated. Ions are produced as a result of the oscillating electrical field between the electrodes and the ion current is subsequently detected. The particular design of the plasma source (Figure 1) avoids excessive dilution of the sample.

An important feature of the BID is the low plasma temperature that does not lead to a significant temperature change (more than  $5 \,^{\circ}$ C) in the quartz tube. This prevents sputtering or material loss from the discharge tube and allows its stable operation for an extended




Figure 1. BID cross-sectional diagram and gas flow. Redrawn from Shinada et al. (2012).




Figure 2. Ionization of the analyte (orange bowl) in a dielectric barrier-discharge sustained in He as plasma gas.

period of time. Furthermore, the low thermal load of the plasma decreases the noise at the detector compared to other high energy plasma sources (Shinada et al. 2012).

In the explanation of the working principles of the GC-barrier discharge ionization detector, the dominant role of the vacuum ultraviolet radiation in the ionization of the sample has to be emphasized. More specifically, it is suggested that the vacuum ultraviolet radiation ( $\Delta \varepsilon = 17.7 \text{ eV}$ ), created by the dielectric barrier discharge and by the excited helium molecules according to the reaction (Figure 2):

$$\text{He}_2(\text{A}^1\Sigma u^+) \rightarrow 2 \text{ He} (1^1S_0)$$

can photoionize the sample (Shinada et al. 2012). Assuming that the ionization of the BID is a photoionization process rather than the combined atomization and ionization that is prevalent in the FID, it is expected that the BID response is concentration-dependent rather than mass flow dependent as with the FID detector.

While it is claimed that the GC-BID system can be applied to a wide range of analytes, there are only few studies to verify its performance. In most cases, the BID has been used to determine a limited number of analytes within a given sample, mostly of low molecular weight and partially oxidized which is the range of analytes for which the BID is said to surpass the FID in its sensitivity and applicability (e.g., measurement of ethanol and water: Weatherly, Woods, and Armstrong 2014; measurement of  $CO_2$  after photocatalytic conversion: Iguchi et al. 2015; a 10 compound set of alkanes, benzenes, alcohols, and FAME after comprehensive two-dimensional GC: Franchina et al. 2015; gaseous formic and acetic acid after needle trap extraction: Ueta et al. 2017, 2018), and greenhouse gases  $CO_2$  and  $N_2O$ : Pascale et al. 2017; as well as permanent gases,  $O_2$  and  $N_2$ : Jo and Kim 2017).

The scope of this study is to systematically evaluate the performance of the commercial GC-BID system by determining a set of standards from various compound classes, such as alcohols, anilines, halogenated hydrocarbons, esters, alkanes, aldehydes, polycyclic aromatic hydrocarbons (PAHs), cyclic compounds, and comparing the results with the well-established flame ionization detector for GC. The analytical figures of merit of the BID are determined and provide a more consistent basis for the assessment of the merits of the dielectric barrier discharge as a GC detector.

#### **Materials and methods**

#### Reagents

All reagents used had a purity of at least 95%. Benzene, butylbenzene, cyclododecane, decane, ethyl acetate, ethylbenzene, isovaleraldehyde, pentachlorophenol, pentylbenzene, 2-propanol, fluoranthene, propylbenzene, *n*-octane, 1-decanol, 1-heptanol, *n*-heptane, cyclododecanone, chrysene, *n*-tetradecane, naphthalene, fluorene, phenanthrene, and pyrene were purchased from Fluka (Seelze, Germany). 1,4-Butanediol, 2-butoxyethyl acetate, 4-chlorophenol, 2,4-dichlorophenol, N,N-diethylamine, N,N-dimethylamine, isopropyl ether, phenol, chlorodecane, chlorododecane, cycloheptane, and cyclodecanone were obtained from Aldrich (Steinheim, Germany). Butylacetate, cyclohexane, cyclodecane, cyclopentane, dodecane, isopropyl acetate, hexane, toluene, 1-pentanol, and cyclopentane were purchased from Sigma-Aldrich (Vienna, Austria). Diethyl ether was purchased from Riedel-de Haën. Benzaldehyde, 1-bromodecane, 1-bromohexane, 1-bromo-3-chloropropane, methanol, 2-methoxy-1-propyl acetate, 1-nonanol, and chloroheptane were obtained from Merck (Vienna, Austria). 1-hexanol, cyclopentanone, and cyclohexanone were obtained from Loba (Fischamend, Austria). Chlorooctane and acenaphthene were from Hüls (Marl, Germany) and Acros Organics (Geel, Belgium), respectively.

All of the reagents were used for the standard solutions. Helium, used as the carrier and discharge gas, had a purity of 99.999% and was purchased by Messer (Austria).

#### Instrumentation

Measurements were performed using a GC-2010 Plus with a BID detector (Tracera) from Shimadzu and a GC-2010 Plus system with flame ionization detector and an AOC-20i autosampler (both Shimadzu, Kyoto, Japan). An RTX-5MS (5% diphenyl/95% dimethyl polysiloxane) GC column of  $30 \text{ m} \times 0.25 \text{ mm} \times 0.25 \mu \text{m}$  from Restek

2826 🛞 M. ANTONIADOU ET AL.

(Bellefonte, PA) and a DB-5MS (5%-diphenyl/95% methylpolysiloxane)  $30 \text{ m} \times 0.25 \text{ mm} \times 0.25 \mu \text{m}$  from Agilent J & W (Middleburg, Netherlands) were employed in this study.

### **Standard solutions**

Standards were prepared in different solvents in order to ensure adequate solubility. Nine groups of standards were prepared at five concentration levels  $(100 \text{ pg}/\mu\text{L}, 1 \text{ ng}/\mu\text{L}, 10 \text{ ng}/\mu\text{L}, 10 \text{ ng}/\mu\text{L}, 1 \mu\text{g}/\mu\text{L})$ . The standards were grouped according to the chemical nature of compounds, namely the anilines, alcohols, aldehydes, halogenated hydrocarbons, phenols, alkanes, cycloalkanes, aromatic compounds, and esters. Alkanes, cyclic and aromatic compounds were dissolved in cyclopentane; the esters in diethyl ether, and the remaining standards in methanol.

Additional standards were analyzed to further test the instrument response. These were the cyclic ketones and cycloalkanes in cyclopentane, chlorinated benzenes, chlorinated alkanes and alcohols in methanol, benzenes, alkanes, and PAHs in dichloromethane. These solutions were prepared at two concentration levels (1 and 100 ng/ $\mu$ L) and an alkane standard (C10-C17) in heptane was prepared at 1 ng/ $\mu$ L.

#### Analytical procedure

Injections of 1  $\mu$ L of each standard were made into the GC-BID and GC-FID system in split mode with split ratio of 20:1. In the optimized method, the BID was operated at a temperature of 300 °C; the FID at a temperature of 250 °C; the injector temperature was 250 °C in both cases; a column flow of 1.08 mL/min, resulting in a carrier gas linear velocity of 25.5 cm/s; purge flow of 3 mL/min; total flow of 25.8 mL/min; inlet pressure of 97.4 kPa; equilibration time of 1 min; and a discharge gas flow of 50 mL/min. The temperature program for the oven was 3 min at 50 °C and increased at 20 °C/min up to 250 °C with a hold time of 4 min. The total analysis time was 17 min.

The second set of standards was analyzed with a similar method but with an initial hold time of only 1 min instead of 3 min which provided an analysis time of 15 min, and with the DB-5MS column. For the PAHs and the alkane standard (C10–C17), the linear ramp of the oven program was extended to  $280 \,^{\circ}$ C and the FID temperature was  $300 \,^{\circ}$ C.


#### **Results and discussion**

#### **Optimization of GC-BID operating conditions**

#### Effect of the septum purge flow rate

The septum purge flow rate is a parameter whose importance is largely underestimated in gas chromatographic analysis. Its optimization is important as a correctly set septum purge flow prevents the contamination of the carrier gas inlet line and removes septum bleed volatiles from the system. On the other hand, an excessive septum purge flow may lead to unintentional dilution of the sample, and hence reduced sensitivity. For this reason, different purge flow rates were examined during analysis. The flow rates ranged from 0.5 to 10 mL/min.

Figure 3 shows the response of a standard with four aromatic compounds. The optimum purge flow rate was 3 mL/min because the lower and higher rates were unable to



**Figure 3.** Effect of the purge flow for aromatic compounds at flow rates of 0, 1, 3, 5, and 10 mL/min. The error bars represent  $\pm$  one standard deviation.


provide greater peak signals—on the low side apparently due to less efficient sample transfer and on the high side due to increased dilution. The current result is in good agreement with the purge flow rate recommended by the manufacturer.

### Effect of the discharge gas flow rate

In order to understand the impact of the system flow rates on the instrument response, we examined the effect of the flow rate of the discharge gas. The instrument can function within a specific flow rate range from relatively low values at which plasma stability may be affected up to ca. 100 mL/min. According to the information provided by the manufacturer, the column flow and discharge gas flow must be adjusted in mutual correspondence: Depending on the particular discharge gas flow rate used, the column flow should be adjusted accordingly. The reason for that apparently lies in the increased pressure at the GC column end when a higher flow rate and thus higher pressure is applied to the detector.

In this work, four flow rates (25, 50, 75, and 100 mL/min) were used for the determination of several compounds at different concentrations and with similar column flow. The evaluation of a set of selected standards at different detector flow rates in Figure 4 shows a noteworthy influence of the detector flow rate on signal intensity. Higher gas flow rate decreases the signal for the studied compounds. Although the discharge gas flow rate of 25 mL/min gave the best response, it was not preferred because it is at the lower end of the acceptable operating range of the instrument, and according to the manufacturer could affect both the reproducibility of measurement as well as have an adverse effect on the instrument long-term stability. Hence, a flow rate of 50 ml/min was chosen for all further measurements, representing a compromise between the anticipated signal stability and observed signal loss. This flow rate is in agreement with the manufacturer's recommended flow rate.

The results from the purge and discharge gas flow rate investigation indicate the high dependence of the instrument response and sensitivity on the flow rate conditions and attest to the importance of optimizing these flow rates.




**Figure 4.** Diagram demonstrating the effect of the discharge gas flow rate for anilines at concentrations of  $10 \text{ ng}/\mu\text{L}$ ,  $100 \text{ ng}/\mu\text{L}$ , and  $1 \mu\text{g}/\mu\text{L}$ ) and gas flow rates of 25, 50, 75, and 100 mL/min. The error bars represent ± one standard deviation. For the higher concentration, the left axis applies, and for the two lower concentrations, the right axis applies.

### Comparison with the FID detector

The measurements reported subsequently are performed under optimum conditions for GC-BID and GC-FID regarding the temperatures and detector flow rates. The same split ratio and temperature program was used for the measurements on both instruments as described in the experimental section. All standards were measured in three replicates on the same day for both detectors.

The chromatograms obtained for both detectors and the investigated sets of standards were perfectly comparable (Figure 5). Notably, no excessive detector-related peak broadening or peak tailing was observed with small exceptions in the case of amine compounds, easily explainable with their high tendency to interact with imperfectly inert surfaces. The retention times for the compounds were consistent between the techniques.

As a means of assessing the performance of the detector under the given chromatographic conditions, the resolution was calculated for the following compounds and their respective neighbor peaks: 1-butanol and isoamyl alcohol; 1-bromo-3-chloropropane and 1-bromohexane; 2,4-dichlorophenol and 4-chlorophenol; butylbenzene and pentylbenzene; butylacetate and 2-methoxy-1-propyl acetate. The resolution factors were in the range from 1.4 to 10.6 for the BID and 2.1-14.1 for the FID and are listed in Table 1. It should be noted that even in the case of the most critical separation, the resolution of the 2,4-dichlorophenol and 4-chlorophenol peaks was still satisfactory. The detector does not



**Figure 5.** Comparison of chromatograms for four groups with the BID and the FID at 100 ng/ $\mu$ L: esters, upper left; alkanes, lower left; cyclic, upper right; and aromatic, lower right. Analytical conditions: injector temperature 250 °C; column flow 1.08 mL/min, 20:1 split; discharge gas flow 50 mL/min; oven initially at 50 °C (3 min) and increased at 20 °C/min to 250 °C (4 min).

Table 1. Chromatographic resolution obtained for the BID and FID using identical separation conditions
for adjacent peak pairs.

Peak pairs	Resolution BID	Resolution FID
1-Butanol/isoamyl alcohol	7.9	8.8
1-Bromo-3-chloropropane/1-bromohexane	7.2	11.4
2,4-Dichlorophenol/4-chlorophenol	1.4	2.1
Butylbenzene/pentylbenzene	10.6	14.1
Butyl acetate/2-methoxy-1-propyl acetate	6.7	9.2

suffer from the introduction of a high concentration of matrix as long as it is sufficiently volatile to not cause further chromatographic problems (e.g., sample carry-over, increased background, peak tailing, or retention time shifts).

Analytical figures of merit were evaluated for all analytes from the corresponding calibration graphs using five concentration levels and are summarized in Table 2. The calibration curves were constructed by plotting the peak area of the signal (y axis) as function of the concentration in the x axis. However, in some cases, the standards could not be detected at the lowest level of the calibration in which case the calibration and calculations were performed with the remaining levels. The coefficients of determination ( $R^2$ ) (Table 2) were typically more than 0.999 for both detectors. The  $R^2$  reached

# 2830 🛞 M. ANTONIADOU ET AL.

Table 2. Comparison of the limit of detection (LOD), limit of quantification (LOQ), and relative stand-
ard deviation (RSD) of the GC-BID and GC-FID.

Analyte   Detector   (pg/s)   (pg/			3σ/S M	Nethod	S/N Method				
NN-Dimethylaniline   BID   315   1041   221   729   11-24 $y = 501e + 65 x - 4499$ 0.9999     NN-Diethylaniline   BID   627   2071   167   553   0.6-4.1 $y = 1.32e + 06 x - 753$ 0.9999     Isovaleraldehyde   BID   197   651   563   1857   0.9-7.3 $y = 1.32e + 06 x - 6751$ 0.9999     Benzaldehyde   BID   202   668   413   1362   0.4-71 $y = 3.80e + 06 x - 65487$ 0.9997     HD   106   351   1184   3006   1.0-21 $y = 1.48798 x - 19452$ 0.9999     HD   106   351   1184   3006   1.0-21 $y = 3.30e + 06 x - 4894$ 0.9999     HD   1059   3493   301   920   2.18 $y = 3.32e + 06 x + 4894$ 0.9999     HD   1059   3493   301   920   4.45.8 $y = 3.32e + 06 x + 4894$ 0.9999     HD   105   3493   301   920   4.45.8 $y = 3.32e + 06 x + 4854$				-		LOQ	RSD (%)	5	
FID1083565641862 $0.9-4.7$ $y = 1.23e + 06 x - 7343$ $0.9995$ NN-DiethylanilineFID1976515631857 $0.9-7.3$ $y = 1.33e + 06 x - 6751$ $0.9999$ IsovaleraldehydeFID2026684131362 $0.6-11$ $y = 3.98e + 06 x - 6751$ $0.9999$ BenzaldehydeFID2026684131362 $0.6-11$ $y = 3.98e + 06 x - 65487$ $0.9999$ BenzaldehydeFID2026684131362 $0.6-11$ $y = 3.98e + 06 x + 4894$ $0.9999$ I-Bromo-3-FID359118615845228 $0.2-13$ $y = 3.902e + 06 x + 4284$ $0.9999$ I-BromohexaneFID34711441111365 $0.2-14$ $y = 3.96e + 05 x + 4284$ $0.9999$ I-BromohexaneFID3009887262396 $0.6 - 1.18$ $y = 3.96e + 05 x - 1134$ $0.9999$ PhenolFID1309887262396 $0.6 - 1.18$ $y = 3.96e + 05 x - 1134$ $0.9999$ ChilorophenolFID19062794663121 $0.4-13$ $y = 77829 x - 10.443$ $0.9999$ ChilorophenolFID1927636014744863 $0.4-72$ $y = 621277 x - 7356$ $0.9999$ ChilorophenolFID1927636014744863 $0.4-72$ $y = 3928 + 06 x - 4530$ $0.9999$ ChilorophenolFID192663666212048 $0.5-53$ $y = 3.728 + 05 x - $	Analyte	Detector	(pg/s)	(pg/s)	(pg/s)	(pg/s)	range	equation	determination R <sup>2</sup>
NN-Dicethylaniline   BID   627   2071   167   553   0.6 - 11   y = 5.45e + 06 x - 1361   0.9999     Isovaleraldehyde   BID   400   1518   140   422   0.6 - 13   y = 3.58e + 06 x + 15720   0.9999     Benzaldehyde   BID   312   1029   420   1386   0.4 - 21   y = 3.58e + 06 x + 15720   0.9999     Benzaldehyde   BID   312   1029   420   1386   0.4 - 21   y = 3.80e + 06 x - 65487   0.9997     1-Bromo-3- chicorpropane   BID   774   2272   89   0.6 - 10   y = 3.920e + 06 x + 2.850   0.9999     1-Bromohexane   BID   774   2533   299   0.6 - 13   y = 3.920e + 06 x + 2.850   0.9999     1-Bromodecane   BID   1059   3493   301   920   0.4 - 31   y = 7.2644 for x + 3810   0.9999     1-Bromodecane   BID   102   217   1111   3050   0.9 - 7.444 for x + 4816   0.9999     2.4-Dichlorophenol   BID   247   1010	N,N-Dimethylaniline								
FID1976515631857 $0.9-7.3$ $y = 1.32 + 0.6 x - 6751$ $0.9996$ IsovaleraldehydeFID2026684131362 $0.6-11$ $y = 3.362 + 0.6 x + 6547$ $0.9997$ BenzaldehydeFID1063511184306 $0.4-21$ $y = 3.980 + 0.6 x - 6547$ $0.9997$ I-Bromo-3-FID1063511184306 $0.4-21$ $y = 3.970 + 0.6 x + 4394$ $0.9999$ I-Bromo-RameFID359118615845228 $0.2-13$ $y = 3.902 + 0.6 x + 4284$ $0.9999$ I-BromohexaneFID347114411113655 $0.2-14$ $y = 3.567.78 x - 3891$ $0.9999$ I-BromohexaneFID3009887262396 $0.9-74$ $y = 7.782.95 x - 10.443$ $0.9999$ I-BromohexaneFID1006279463121 $0.4-13$ $y = 3.96e + 0.6 x - 1344$ $0.9999$ I-BromohexaneFID1006279463121 $0.4-72$ $y = 64.298 x - 67.59$ $0.9999$ I-DromohexaneFID102636014744863 $0.4-72$ $y = 64.298 x - 67.59$ $0.9999$ I-DromohexaneFID1027636014744863 $0.4-72$ $y = 3.12e + 0.6 x - 1344$ $0.9999$ I-DromohexaneFID110223776310.5-19 $y = 3.12e + 0.6 x - 1344$ $0.9999$ I-DromohexaneFID11022361201234 $0.6-5.3$ $y = 3.98e + 0.6 x - 1344$								-	
isovaleraldehyde FID   400   4518   140   462   0.6-11   y = 3,562 + 0.6 x + 15730   0.9999     Benzaldehyde FID   010   312   1029   420   1386   0.4-21   y = 3,802 + 0.6 x - 65437   0.997     1-Bromo-3- chicorpropane   BID   22   272   288   0.6-10   y = 2,972 + 0.6 x + 4984   0.9999     1-Bromohexane   BID   774   2533   299   987   0.6-13   y = 3,022 + 0.6 x + 2,850   0.9999     1-Bromohexane   BID   774   2533   299   0.974   y = 3,022 + 0.6 x + 2,850   0.9993     1-Bromodecane   BID   1059   3493   301   92   0.4-5.8   y = 3,462 + 0.6 x + 1314   0.9999     1-Bromodecane   BID   427   1410   351   1159   0.3-5.9   y = 3,162 + 0.6 x - 1344   0.9999     2.4-Dichlorophenol   BID   342   104   351   1159   0.3-5.9   y = 3,162 + 0.6 x - 1344   0.9999     2.4-Dichlorophenol   BID   420   1411 <th>N,N-Diethylaniline</th> <th></th> <th></th> <th></th> <th></th> <th></th> <th></th> <th></th> <th></th>	N,N-Diethylaniline								
FID20266841313620.4-38 $y = 494098 x - 56347$ 0.9999BenzaldehydeBID31210/210/80.6-10 $y = 3.06 + 6.06 \times - 55487$ 0.99971-Brono-3-BID6.2920752728980.6-10 $y = 2.97e + 0.6 \times + 4984$ 0.99991-Brono-A-BID77425532999870.6-35 $y = 3.02e + 0.6 \times + 4895$ 0.99991-BromohexaneBID77425532999870.6-35 $y = 3.56778 \times - 3891$ 0.99991-BromodecaneBID105934933019920.4-5.8 $y = 3.56e + 0.6 \times + 11314$ 0.99991-BromodecaneBID105934933019920.4-5.4 $y = 3.36e + 0.6 \times - 13484$ 0.9999PhenolBID186614281920.4-7.4 $y = 3.56e + 0.6 \times - 13444$ 0.99992-4'DichlorophenolBID102744015111590.3-5.9 $y = 3.14e + 0.6 \times - 13944$ 0.99992-4'DichlorophenolBID1227451304281.4-4.7 $y = 2.42e + 0.6 \times - 1394$ 0.99992-4'DichlorophenolBID2267451304281.4-4.7 $y = 3.37e + 0.6 \times - 1394$ 0.99992-A'DichlorophenolBID2267451304281.4-4.7 $y = 2.42e + 0.6 \times - 82328$ 0.99992-A'DichlorophenolBID2267451304281.4-4.7 $y = 3.37e + 0.6 \times - 32828$ 0.99992-Bro	Is a second s								
Benzaldehyde IDBID31210294201360.4-21 $y = 3.80e + 06 \times - 65.470.9971-Bromo-3-chloropropaneBID62920752728980.6-10y = 2.97e + 06 \times + 49840.99991-Bromo-A-chloropropaneFID3591186158452280.2-13y = 3.95005 \times - 2.9870.99931-BromohexaneBID77425532999870.6-35y = 3.02e + 06 \times + 43840.99991-BromodecaneBID147114411113650.2-14y = 5.5678 \times - 38910.99941-BromodecaneBID15034933019920.4-5.8y = 3.36e + 06 \times + 131440.999924-DichlorophenolBID1866142819260.1-1.8y = 3.56e + 06 \times - 141440.99992.4-DichlorophenolBID13411022317610.5-19y = 3.12e + 06 \times + 15340.99992.4-DichlorophenolBID2267463121.440.45 - 123780.99992-PropanolBID126630144647710.7-42y = 6321277 \times -73560.99992-BropanolBID1267463120.49y = 3.37e + 06 \times + 5300.99991-ButanolBID12860014461710.7-42y = 632127 \times -73560.99991-ButanolBID1287206441260.4-72y = 3.37e + 06 \times + 53030.99991-ButanolB$	isovaleraidenyde								
Ib106351118439061.0-21 $y = 1.43798 \times - 194220.9991-Bromo-SaneFID3591186158452280.2-13y = 3.9050 \times - 29870.99991-BromohexaneFID3471144111136550.2-13y = 3.9050 \times - 29870.99991-BromohexaneFID3471144111136550.2-14y = 3.56778 \times - 38910.99991-BromodecaneBID105934933019920.4-5.8y = 3.56778 \times - 38910.99991-BromodecaneBID105934933019920.4-5.8y = 3.56+766 \times - 141840.99991-BromohexaneBID126764621210.4-31y = 7.78259 \times - 104430.99992.4-DichlorophenolBID1276360147448630.4-72y = 3.12e+766 \times - 5300.99992.4-DichlorophenolBID33411022317610.5-19y = 3.17e+766 \times + 5300.99992-PropanolBID2267451304281.4-4.7y = 3.32e+766 \times - 32020.99991-ButanolBID22870764422890.5-65y = 3.39e+76 \times - 33020.99991-ButanolBID54816072237380.8-11y = 3.13e+76 \times + 1.5040.99991-ButanolBID54816072237380.8-11y = 3.13e+76 \times + 1.5040.99991-ButanolBID5481607223$	Ponzaldohudo							-	
1-Bromo-3- chloropropane   BD   629   2075   272   898   0.6-10 $y = 2.97e + 06 x + 4984$ 0.9999     1-Bromohexane   BID   774   2553   299   987   0.6-35 $y = 3.02e + 06 x + 4384$ 0.9999     1-Bromohexane   BID   774   2553   299   987   0.6-35 $y = 3.02e + 06 x + 1381$ 0.9994     1-Bromodecane   BID   1059   3493   301   992   0.4-5.8 $y = 3.32e + 06 x + 1314$ 0.9999     Phenol   BID   1057   746   312   0.4-31   0.9999 $y = 7326446 x - 4816$ 0.9999     2.4-Dichlorophenol   BID   427   1410   351   1159   0.3-5.9 $y = 3.12e + 06 x - 14134$ 0.9999     4-Chlorophenol   BID   334   1102   231   761   0.5-19 $y = 3.12e + 06 x - 15236$ 0.9999     2-Propanol   BID   334   1102   231   741   0.7-42 $y = 4.0225 x - 6579$ 0.9990     2-Bropanol   BID	benzaluenyue								
chloropropane   FID   359   1186   1524   5228   0.2-13   y = 395005 x - 2987   0.9993     1-Bromodecane   BID   774   2553   299   987   0.6-35   y = 302e + 06 x + 2850   0.9999     1-Bromodecane   BID   1059   3493   301   92   0.4-5.8   y = 3362 + 06 x + 11314   0.9999     Phenol   BID   1059   3493   301   92   0.4-5.8   y = 3782 + 06 x + 11314   0.9999     Phenol   BID   136   614   281   92   0.4-131   y = 778259 x - 10443   0.9999     2.4-Dichlorophenol   BID   334   1102   231   761   0.5-19   y = 3172 + 06 x + 2530   0.9999     2.4-Chorophenol   BID   226   745   130   428   1.4-4.7   y = 3428 + 06 x - 21324   0.9999     2.4-Chorophenol   BID   226   745   130   428   1.6-5.13   y = 31521 x - 8632   0.9999     2.4-Dichlorophenol   BID   248	1-Bromo-3-							,	
FID359118615845228 $0.2-13$ $y = 395005 x - 2987$ $0.9993$ 1-BromohexaneBID7742553299 $0.6-35$ $y = 356778 x - 3891$ $0.9994$ 1-BromohexaneBID1059343301992 $0.4-58$ $y = 3.48e + 06 x + 11314$ $0.9999$ PhenolBID136614281926 $0.9-74$ $x = 72464 x - 4816$ $0.9993$ PhenolBID136614281926 $0.1-18$ $y = 73269 x - 10443$ $0.9999$ 2.4-DichlorophenolBID3271303511159 $0.3-59$ $y = 3.14e + 06 x - 1334$ $0.9999$ 4-ChlorophenolBID3341102231761 $0.5-19$ $y = 3.12e + 06 x + 5636$ $0.9999$ 2-PropanolBID2267451304281.4-4.7 $y = 3.42e + 06 x - 82938$ $0.999$ 1-ButanolBID528720644288 $0.5-65$ $y = 3.5131 x - 8632$ $0.999$ 1-ButanolBID5481807223738 $0.8-11$ $y = 3.12e + 06 x - 3202$ $0.9999$ 1-ButanolBID5481807223718 $0.8-11$ $y = 3.13e + 06 x - 24660.99991-ButanolBID1464721271.1-2y = 3.5675 x + 5880.99991-ButanolBID54818072237180.9-14y = 3.13e + 06 x - 32020.99991-ButanolBID146477157<$			02)	2075	272	070	0.0 10	y = 2.57C + 00 x + 4504	0.9999
1-BromohexaneBID77425532999870.6-35 $y' = 3.02 + 06 x + 2850$ 0.99991-BromodecaneBID10593433019920.4-5.8 $y = 3.48 + 106 x + 11314$ 0.99931-BromodecaneBID105934333019920.4-5.8 $y = 3.48 + 106 x + 11314$ 0.9999PhenolBID10662794631210.4-31 $y = 778259 x - 10443$ 0.99992.4-DichlorophenolBID427141035111590.3-59 $y = 3.14e + 06 x - 1394$ 0.99994-ChicophenolBID34311022317110.7-42 $y = 621277 x - 7356$ 0.99934-ChicophenolBID2267451304281.4-4.7 $y = 3.42e + 06 x - 2303$ 0.9992-PropanolBID2267451304281.4-4.7 $y = 3.42e + 06 x - 3202$ 0.99992-PropanolBID2267451304281.4-4.7 $y = 3.42e + 06 x - 3202$ 0.99991-ButanolBID42914161384551.0-53 $y = 2.98e + 06 x - 3202$ 0.99991-ButanolBID54818072233600.5-65 $y = 3.47e + 06 x + 1504$ 0.99991-A-ButanediolBID1073333651200.5-14 $y = 3.57e + 588$ 0.99991-A-ButanediolBID1073331200.5-13 $y = 3.49e + 06 x - 1504$ 0.99991-A-ButanediolBID1076		FID	359	1186	1584	5228	0.2-13	$v = 395005 \ x - 2987$	0.9993
FID347114411113665 $0.2-14$ $y = 556778 \times -38911$ 0.99941-BromodecaneBID10593493301920 $0.4-5.8$ $y = 3.86 + 06 \times +11314$ 0.9999PhenolBID186614281926 $0.1-1.8$ $y = 3.36e + 06 \times +11814$ 0.99992.4-DichlorophenolBID42714103511159 $0.3-5.9$ $y = 3.14e + 06 \times -1394$ 0.99992.4-DichlorophenolBID3341102231761 $0.5-19$ $y = 3.17e + 06 \times + 5630$ 0.99992.4-DichlorophenolBID3341102231761 $0.5-19$ $y = 3.17e + 06 \times + 5630$ 0.99992.4-PropanolBID2267451304281.4-4.7 $y = 3.42e + 06 \times - 3293$ 0.9992.PropanolBID2267451304281.4-4.7 $y = 3.42e + 06 \times - 3293$ 0.9991-8utanolBID42914161384551.0-5.3 $y = 3.5131 \times - 8632$ 0.9991-somyl alcoholBID173336512030.9-1.4 $y = 3.13e + 06 \times - 3466$ 0.99991-6utanolBID10735336512030.9-1.4 $y = 3.19e + 06 \times - 3456$ 0.99991-6utanolBID176487222113121.1-1 $y = 3.49e + 06 \times - 3456$ 0.99991-6utanolBID176487225113200.4-7.2 $y = 3.19e + 06 \times - 34576$ 0.99991-6utanolBID147	1-Bromohexane				299				
FID3009887262396 $0.9-7.4$ $y = 7.26446 x - 4816$ $0.9993$ PhenolBID186614281926 $0.1-1.8$ $y = 3.96 + 06 x - 14184$ $0.9990$ 2.4-DichlorophenolBID42714103511159 $0.3-5.9$ $y = 3.14e + 06 x - 1394$ $0.9999$ 4-ChlorophenolBID3341102231761 $0.5-19$ $y = 3.17e + 06 x + 5630$ $0.9999$ 4-ChlorophenolBID226745130428 $1.4-4.7$ $y = 3.42e + 06 x - 8739$ $0.9999$ 2-PropanolBID226745130428 $1.4-4.7$ $y = 3.42e + 06 x - 8238$ $0.9999$ 1-ButanolBID4291416138455 $1.0-5.3$ $y = 2.98e + 06 x - 3202$ $0.9999$ 1-ButanolBID5481807223738 $0.8-11$ $y = 3.13e + 06 x - 2486$ $0.9999$ 1-butanolBID5481807223738 $0.8-11$ $y = 3.13e + 06 x - 2486$ $0.9999$ 1-AbutanediolBID1733533651203 $0.9-1.4$ $y = 3.15e + 06 x + 1.504$ $0.9999$ 1-ctanolBID14764872251252 $0.5-6.5$ $y = 4.76103 x - 2.6876$ $0.9995$ 1-ctanolBID141313882272 $0.4-2.0$ $y = 7.09e + 06 x - 4876$ $0.9999$ 1-dctanolBID14764872251270 $0.3-6.0$ $y = 1.23e + 06 x - 4579$ $0.9999$ <t< th=""><th></th><th>FID</th><th>347</th><th>1144</th><th>1111</th><th>3665</th><th></th><th></th><th></th></t<>		FID	347	1144	1111	3665			
Phenol   BID   B86   614   281   926   0.1-1.8   y = 3.96e + 05.x - 14184   0.9999     2.4-Dichlorophenol   BID   427   1410   311   10.4-31   y = 778259 x - 10443   0.9999     2.4-Dichlorophenol   BID   127   1410   311   159   0.3-5.9   y = 3.126 + 05.x - 1394   0.9999     4-Chlorophenol   BID   334   1102   231   761   0.5-19   y = 3.17e + 05.x + 530   0.9999     2-Propanol   BID   226   745   130   428   1.4-4.7   y = 3.42e + 06.x - 82938   0.9999     1-Butanol   BID   226   745   130   428   1.4-4.7   y = 3.42e + 06.x - 82938   0.9999     1-statanol   BID   226   745   130   428   1.4-4.7   y = 3.5131 x = 6362   0.9999     1-statanol   BID   235   1073   644   2126   0.4-7.2   y = 3.928 + 05.x + 1.504   0.9999     1-statanol   BID   1025   3355	1-Bromodecane	BID	1059	3493	301	992	0.4-5.8	y = 3.48e + 06 x + 11314	0.9999
FID19062794631210.4–31 $y = 772298 + 10443$ 0.99992.4-DichlorophenolBID427141035111590.3–59 $y = 3.14e + 06x - 1394$ 0.99994-ChlorophenolBID33411022317610.5–19 $y = 3.17e + 06x + 5630$ 0.99994-ChlorophenolBID2267451304281.4–4.7 $y = 3.42e + 06x - 82938$ 0.99992-PropanolBID2267451304281.4–4.7 $y = 3.42e + 06x - 82938$ 0.9991-ButanolBID2267451304281.4–4.7 $y = 3.42e + 06x - 82938$ 0.99991-ButanolBID2267451304281.4–4.7 $y = 3.42e + 06x - 82938$ 0.99991-SeutonolBID12872069422890.5–6.5 $y = 476348x - 3727$ 0.9994Isoamyl alcoholBID54818072237380.8–11 $y = 3.13e + 06x - 1304$ 0.99991,4-ButanediolBID10735336512030.9–14 $y = 3.15e + 106x + 1.504$ 0.99951,0-CtanolBID147648722518290.6–53 $y = 3.49e + 06x + 1504$ 0.99951-OctanolBID74744871257131.1–12 $y = 404225x - 3513$ 0.9996HeptaneBID74113882720.4–20 $y = 7.09e + 06x + 45876$ 0.9995DoceaneBID7542487125 <th></th> <th>FID</th> <th>300</th> <th>988</th> <th>726</th> <th>2396</th> <th>0.9–7.4</th> <th><i>y</i> = 726446 <i>x</i> - 4816</th> <th>0.9993</th>		FID	300	988	726	2396	0.9–7.4	<i>y</i> = 726446 <i>x</i> - 4816	0.9993
2.4-DichlorophenolBID427141035111590.3-5.9 $y = 3.14 \pm 0.6 x - 1394$ 0.99994-ChlorophenolBID33411022317610.5-19 $y = 3.17e \pm 0.6 x - 1394$ 0.99992-PropanolBID2267451304281.4-7 $y = 621277 x - 7356$ 0.99992-PropanolBID2267451304281.4-7 $y = 3.42e \pm 0.6 x - 82938$ 0.9991-ButanolBID42014161384551.0-5.3 $y = 2.98e \pm 0.6 x - 3202$ 0.99991-ButanolBID54818072237380.8-11 $y = 3.13e \pm 0.6 x - 2486$ 0.99991-ButanolBID54818072237380.8-11 $y = 3.13e \pm 0.6 x - 2486$ 0.99991,4-ButanediolBID54818072237380.8-11 $y = 3.13e \pm 0.6 x - 2486$ 0.99991,4-ButanediolBID117335335512030.9-14 $y = 3.15e \pm 0.6 x + 1.504$ 0.99981,4-ButanediolBID147648722518290.6-53 $y = 3.49e \pm 0.6 x + 1504$ 0.99991-OctanolBID147648722518290.6-53 $y = 3.49e \pm 0.6 x + 45876$ 0.99991-OctanolBID75424871254130.1-13 $y = 7.55e \pm 0.6 x - 45639$ 0.9996DecaneBID75424871254130.1-13 $y = 7.55e \pm 0.6 x - 45639$ 0.9993Decane<	Phenol		186	614	281	926	0.1–1.8	y = 3.96e + 06 x - 14184	0.9999
FID19276360147448630.4-72 $y = 62127x - 7356$ 0.99994-ChlorophenolBID33411022317610.5-19 $y = 3.17e + 06x + 5630$ 0.99992-PropanolBID2267451304281.4-4.7 $y = 3.42e + 06x - 82938$ 0.9991-ButanolBID42914161384551.0-5.3 $y = 2.98e + 06x - 82923$ 0.99991-ButanolBID42914161384551.0-5.3 $y = 2.98e + 06x - 3202$ 0.99991-ButanolBID42914161384551.0-5.3 $y = 2.98e + 06x - 3202$ 0.99991-ButanolBID42914161384551.0-5.3 $y = 2.98e + 06x - 3202$ 0.99991-ButanolBID14718072237380.8-11 $y = 3.13e + 06x - 3202$ 0.99991-A-ButanediolBID14748702273780.8-11 $y = 3.13e + 06x - 3202$ 0.99991-A-ButanediolBID14748722518290.6-5.3 $y = 3.49e + 06x + 1.504$ 0.99991-A-ButanediolBID147648722518290.6-5.4 $y = 3.49e + 06x + 1.504$ 0.99991-OctanolBID4111358822720.4-2.0 $y = 7.09e + 06x - 45876$ 0.99951-OctanolBID4111358822720.4-2.0 $y = 7.09e + 06x - 45876$ 0.9999DecaneBID76224871									
4-Chlorophenol   BID   334   1102   231   761   0.5-19   y = 3.17e+06x+5630   0.9999     2-Propanol   BID   226   600   1446   4771   0.7-42   y = 642988 x - 6759   0.9996     2-Propanol   BID   226   745   130   428   1.4-4.7   y = 3.42e+06x - 82938   0.999     1-Butanol   BID   429   1416   138   455   1.0-5.3   y = 2.98e+06x - 3202   0.9999     Isoamyl alcohol   BID   218   720   694   2289   0.5-6.5   y = 3.43e+06x - 2486   0.9999     Isoamyl alcohol   BID   107   353   365   1203   0.9-1.4   y = 3.15e+06x+1.504   0.9999     1.4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e+06x+1.504   0.9999     1.0-Ctanol   BID   1476   4872   227   0.4-2.0   y = 7.09e+06x - 45876   0.9999     1.0-Ctanol   BID   754   2487   125 <t< th=""><th>2.4-Dichlorophenol</th><th></th><th></th><th></th><th></th><th></th><th></th><th>-</th><th></th></t<>	2.4-Dichlorophenol							-	
FID   182   600   1446   471   0.7-42   y = 642988 x - 6759   0.9996     2-Propanol   BID   226   745   130   428   1.4-4.7   y = 3.42e + 06 x - 82938   0.999     1-Butanol   BID   429   1416   138   455   1.0-5.3   y = 2.98e + 06 x - 3202   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.13e + 06 x - 2486   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.13e + 06 x - 1.504   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.13e + 06 x + 1.504   0.9998     1,4-Butanediol   BID   107   353   365   1203   0.9-1.1   y = 3.13e + 06 x + 1.504   0.9998     1,4-Butanediol   BID   1476   4872   251   829   0.6-5.3   y = 3.49e + 06 x + 1.504   0.9998     1,4-Butanediol   BID   2105   9106									
2-Propanol   BID   226   745   130   428   1.4-4.7   y = 3.42e + 06 x - 82938   0.999     HD   5108   16856   621   2048   0.6-5.8   y = 3.5511 x - 8632   0.999     I-Butanol   BID   429   1416   138   455   1.0-5.3   y = 2.98e + 06 x - 3202   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.15e + 06 x - 2.486   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e + 06 x + 1.504   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e + 06 x + 1.504   0.9999     1,0-Ctanol   BID   1476   4872   221   62.72   y = 476103 x - 2655   0.9998     Heptane   BID   2760   9106   344   1134   0.5-4.1   y = 676103 x - 2655   0.9999     Decane   BID   754   2487   125   413 <t< th=""><th>4-Chlorophenol</th><th></th><th></th><th></th><th></th><th></th><th></th><th></th><th></th></t<>	4-Chlorophenol								
FID   5108   16856   621   2048   0.6-5.8   y = 355131 x - 8632   0.999     1-Butanol   BID   429   1416   138   455   1.0-5.3   y = 2.984 + 06 x - 3202   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.13e + 06 x - 2486   0.9999     1.4-Butanediol   BID   107   353   355   1203   0.9-14   y = 3.15e + 06 x + 1.504   0.9999     1.4-Butanediol   BID   107   353   355   1203   0.9-14   y = 3.15e + 06 x + 1504   0.9999     1.4-Butanediol   BID   1476   4872   251   829   0.6-5.3   y = 3.987 x + 588   0.9998     Heptane   BID   411   1358   82   272   0.4-2.0   y = 7.09e + 06 x - 45876   0.9999     Decane   BID   752   413   0.1-1.3   y = 7.55e + 06 x - 45639   0.9996     Cyclohexane   BID   637   2103   79   260   0.4-3.7	2 Duemental								
1-Butanol   BID   429   1416   138   455   1.0-5.3   y = 2.98e + 06 x - 3202   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.13e + 06 x - 2486   0.9999     Isoamyl alcohol   BID   548   1807   223   738   0.8-11   y = 3.13e + 106 x - 2486   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e + 06 x + 1.504   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e + 06 x + 1.504   0.9999     1,4-Butanediol   BID   1476   4872   251   829   0.6-5.3   y = 3.49e + 06 x + 1504   0.9998     Heptane   BID   716   4872   251   709   0.3-6.0   y = 7.09e + 06 x - 45876   0.99993     Decane   BID   754   2487   125   413   0.1-1.3   y = 7.55e + 06 x - 45639   0.9993     Dodecane   BID   562	2-Propanoi							•	
FID21872069422890.5-6.5 $y = 4.76348 x - 3727$ 0.9994Isoamyl alcoholBID54818072237380.8-11 $y = 3.13e + 06 x - 2486$ 0.99991,4-ButanediolBID10735336512030.9-1.4 $y = 3.15e + 06 x + 1.504$ 0.99991,4-ButanediolBID10735336512030.9-1.4 $y = 3.15e + 06 x + 1.504$ 0.99991-OctanolBID147648722518290.6-5.3 $y = 3.49e + 06 x + 1504$ 0.9998HeptaneBID4111358822720.4-2.0 $y = 7.09e + 06 x - 45876$ 0.9995DecaneBID75424871254130.1-1.3 $y = 7.55e + 06 x - 45639$ 0.9996DecaneBID56218561504960.3-13 $y = 7.55e + 06 x - 45639$ 0.9997DodecaneBID6372103792600.4-3.7 $y = 3.99e + 06 x - 12696$ 0.9993CyclohexaneBID6372103792600.4-3.7 $y = 3.99e + 06 x - 16656$ 0.9993CyclodecaneBID6372103792600.4-3.7 $y = 3.39e + 06 x - 2164$ 0.9999FID1504962287530.9-2.3 $y = 6.49850 x - 3228$ 0.9993CyclodecaneBID637210379200.4-3.7 $y = 3.39e + 06 x - 16656$ 0.9993CyclodecaneBID1304291972	1-Rutanol								
Isoamyl alcohol FID   BID FID   548 325   1073 1073   644 2126   0.4-7.2 0.4-7.2   y = 3.15e + 06 x + 1.504 y = 3.15e + 06 x + 1.504   0.9999     1,4-Butanediol BID   107   353   365   120   0.4-7.2   y = 3.15e + 06 x + 1.504   0.9999     1-Octanol FID   198   654   659   2173   1.1-12   y = 404225 x - 3513   0.9995     1-Octanol FID   2760   9106   344   1134   0.5-4.1   y = 7.6103 x - 2655   0.9998     Heptane FID   926   3055   215   709   0.3-6.0   y = 1.23e + 06 x - 1504   0.9999     Decane FID   920   660   260   859   0.4-7.5   y = 7.51e + 06 x - 12046   0.9993     Dodecane FID   78   258   261   860   0.3-17   y = 1.52e + 06 x - 12046   0.9993     Cyclohexane FID   715   496   228   753   0.9-2.3   y = 6.49e + 06 x - 16656   0.9993     Cyclodecane FID   150   496   228   753   0.9-2.3   y = 6.49e + 06 x - 10544								,	
FID   325   1073   644   2126   0.4-7.2   y = 395875 x + 588   0.9999     1,4-Butanediol   BID   107   353   365   1203   0.9-1.4   y = 3.15e + 06 x + 1.504   0.9999     1-Octanol   BID   1476   4872   251   829   0.6-5.3   y = 3.49e + 06 x + 1504   0.9998     Heptane   BID   4476   4872   251   829   0.6-5.3   y = 3.49e + 06 x + 1504   0.9998     Heptane   BID   411   1358   82   272   0.4-2.0   y = 7.09e + 06 x - 45876   0.9994     Decane   BID   754   2487   125   413   0.1-1.3   y = 7.55e + 06 x - 45639   0.9994     Dodecane   BID   562   1856   150   496   0.3-13   y = 7.52e + 06 x - 12046   0.9993     Cyclohexane   BID   657   2188   100   0.3-17   y = 1.69e + 06 x - 12695   0.9994     Cyclohecane   BID   657   2168   129   426   0	Isoamvi alcohol								
1,4-Butanediol BID 107 353 365 1203 0.9-1.4 y = 3.15e + 06 x + 1.504 0.9999   1-Octanol BID 1476 4872 251 829 0.6-5.3 y = 3.49e + 06 x + 1504 0.9998   Heptane BID 411 1358 82 272 0.4-2.0 y = 7.09e + 06 x - 45876 0.9998   Heptane BID 411 1358 82 272 0.4-2.0 y = 7.09e + 06 x - 45876 0.9999   Decane BID 754 2487 125 413 0.1-1.3 y = 7.55e + 06 x - 12046 0.9993   Dodecane BID 562 1856 150 496 0.3-17 y = 1.52e + 06 x - 12046 0.9993   Cyclohexane BID 562 1856 150 496 0.3-17 y = 1.69e + 06 x - 12695 0.9994   Cyclohexane BID 657 2168 129 426 0.3-2.3 y = 6.48950 x - 3228 0.9993   Cyclododecane BID 657 2168 129 426 0.3-2.3 y = 6.32e + 06 x - 10544 0.9999									
FID19865465921731.1-12 $y = 404225 x - 3513$ 0.99951-OctanolBID147648722518290.6-5.3 $y = 3.49e + 06 x + 1504$ 0.9998FID2760910634411340.5-4.1 $y = 676103 x - 2655$ 0.9998HeptaneBID4111358822720.4-2.0 $y = 7.09e + 06 x - 45876$ 0.9995FID92630552157090.3-6.0 $y = 1.23e + 06 x - 9190$ 0.9994DecaneBID75424871254130.1-1.3 $y = 7.55e + 06 x - 12046$ 0.9993DodecaneBID56218561504960.3-13 $y = 7.61e + 06 x - 38652$ 0.9997FID782582618600.3-17 $y = 1.52e + 06 x - 12046$ 0.9993CyclohexaneBID6372103792600.4-3.7 $y = 3.99e + 06 x - 12695$ 0.9994CyclodocaneBID65721681294260.3-2.3 $y = 6.48950 x - 3228$ 0.9993CyclododecaneBID65721681294260.3-2.3 $y = 6.33e + 06 x - 10544$ 0.9999FID1304291053460.2-37 $y = 97133 x - 7485$ 0.9993CyclododecaneBID6332088973200.1-5.0 $y = 4.06e + 06 x - 4874$ 0.9999FID250382611836030.6-3.3 $y = 1.26e + 06 x - 4874$ 0.9999Tolue	1.4-Butanediol							-	
1-Octanol BID 1476 4872 251 829 0.6-5.3 y = 3.49e + 06 x + 1504 0.9998   Heptane BID 2760 9106 344 1134 0.5-4.1 y = 6.76103 x - 2655 0.9998   Heptane BID 411 1358 82 272 0.4-2.0 y = 7.09e + 06 x - 45876 0.9995   Decane BID 754 2487 125 413 0.1-1.3 y = 7.55e + 06 x - 45639 0.9996   Dodecane BID 754 2487 125 413 0.1-1.3 y = 7.55e + 06 x - 45639 0.9997   Cyclohexane BID 562 1856 150 496 0.3-13 y = 7.51e + 06 x - 12695 0.9994   Cyclohexane BID 637 2103 79 260 0.4-3.7 y = 3.99e + 06 x - 16656 0.9993   Cyclodecane BID 637 2103 79 260 0.4-3.7 y = 6.38e + 06 x - 16656 0.9993   Cyclodecane BID 637 2168 129 426 0.3-2.3 y = 6.48950 x - 3228 0.9993	.,								
Heptane   BID   411   1358   82   272   0.4-2.0   y = 7.09e + 06 x - 45876   0.9995     Decane   BID   754   2487   125   709   0.3-6.0   y = 7.55e + 06 x - 45639   0.9996     FID   200   660   260   859   0.4-7.6   y = 1.52e + 06 x - 12046   0.9993     Dodecane   BID   562   1856   150   496   0.3-13   y = 7.61e + 06 x - 13265   0.9993     Cyclohexane   BID   637   2103   79   260   0.4-3.7   y = 3.99e + 06 x - 16656   0.9993     Cyclodecane   BID   637   2168   129   426   0.3-2.3   y = 6.04e + 06 x - 10544   0.9999     Cyclododecane   BID   637   1168   129   426   0.2-34   y = 6.32e + 06 x - 2164   0.9999     Cyclododecane   BID   130   429   105   346   0.2-34   y = 6.32e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2	1-Octanol	BID	1476	4872	251	829	0.6-5.3		
FID   926   3055   215   709   0.3-6.0   y = 1.23e + 06 x - 9190   0.9994     Decane   BID   754   2487   125   413   0.1-1.3   y = 7.55e + 06 x - 45639   0.9996     Dodecane   BID   562   1856   150   496   0.3-13   y = 7.51e + 06 x - 12046   0.9993     Dodecane   BID   562   1856   150   496   0.3-17   y = 1.52e + 06 x - 12046   0.9993     Cyclohexane   BID   637   2103   79   260   0.4-3.7   y = 3.99e + 06 x - 16656   0.9993     Cyclodecane   BID   637   2103   79   260   0.4-3.7   y = 3.99e + 06 x - 16556   0.9993     Cyclodecane   BID   657   2168   129   426   0.3-2.3   y = 6.48950 x - 3228   0.9993     Cyclododecane   BID   130   429   105   346   0.2-3.3   y = 6.33e + 06 x - 164   0.9999     FID   3634   1194   295   972   2.3-6.		FID	2760	9106	344	1134	0.5-4.1	$y = 676103 \ x - 2655$	0.9998
Decane   BID   754   2487   125   413   0.1–1.3   y = 7.55e + 06 x - 45639   0.9996     FID   200   660   260   859   0.4–7.6   y = 1.52e + 06 x - 12046   0.9993     Dodecane   BID   562   1856   150   496   0.3–17   y = 7.61e + 06 x - 12695   0.9994     Cyclohexane   BID   637   2103   79   260   0.4–3.7   y = 3.99e + 06 x - 16656   0.9993     Cyclodecane   BID   637   2163   496   228   753   0.9–2.3   y = 6.48950 x - 3228   0.9993     Cyclodecane   BID   657   2168   129   426   0.3–2.3   y = 6.48950 x - 3228   0.9993     Cyclododecane   BID   475   1567   419   1384   0.2–7.7   y = 971333 x - 7485   0.9999     FID   3634   11994   295   972   2.3–6.2   y = 1.45e + 06 x - 4874   0.9999     FID   3634   1199   295   972   2.3–6.2	Heptane	BID	411	1358	82	272	0.4-2.0	y = 7.09e + 06 x - 45876	0.9995
FID200660260859 $0.4-7.6$ $y = 1.52e + 06 x - 12046$ $0.9993$ DodecaneBID5621856150496 $0.3-13$ $y = 7.61e + 06 x - 38652$ $0.9997$ FID78258261860 $0.3-17$ $y = 1.69e + 06 x - 12695$ $0.9993$ CyclohexaneBID637210379260 $0.4-3.7$ $y = 3.99e + 06 x - 16556$ $0.9993$ CyclodecaneBID6572168129426 $0.3-2.3$ $y = 6.04e + 06 x - 10544$ $0.9999$ CyclododecaneBID47515674191384 $0.2-7.7$ $y = 971333 x - 7485$ $0.9993$ CyclododecaneBID130429105346 $0.2-34$ $y = 6.33e + 06 x - 2164$ $0.9999$ FID363411994295972 $2.3-6.2$ $y = 1.45e + 06 x - 8385$ $0.9994$ BenzeneBID633208897320 $0.1-50$ $y = 4.06e + 06 x - 4874$ $0.9999$ TolueneBID4771574124408 $0.1-10$ $y = 4.40e + 06 x + 1766$ $0.9999$ EthylbenzeneBID6132025140462 $0.4-2.3$ $y = 4.60e + 06 x - 2355$ $0.9999$ ButylbenzeneBID6132025140462 $0.3-3.3$ $y = 5.86e + 06 x + 2456$ $0.9999$ FID1464803691216 $0.9-1.6$ $y = 1.32e + 06 x - 5377$ $0.9999$ ButylbenzeneBID7492472									0.9994
DodecaneBID5621856150496 $0.3-13$ $y = 7.61e + 06x - 38652$ $0.9997$ CyclohexaneBID637210379260 $0.4-3.7$ $y = 1.69e + 06x - 12695$ $0.9994$ CyclohexaneBID637210379260 $0.4-3.7$ $y = 3.99e + 06x - 16656$ $0.9993$ CyclodecaneBID6572168129426 $0.3-2.3$ $y = 6.04e + 06x - 3228$ $0.9993$ CyclododecaneBID47515674191384 $0.2-7.7$ $y = 971333x - 7485$ $0.9993$ CyclododecaneBID130429105346 $0.2-34$ $y = 6.38e + 06x - 2164$ $0.9999$ FID363411994295972 $2.3-6.2$ $y = 1.45e + 06x - 8385$ $0.9994$ BenzeneBID633208897320 $0.1-5.0$ $y = 4.06e + 06x - 492$ $0.9999$ TolueneBID4771574124408 $0.1-10$ $y = 4.40e + 06x + 1766$ $0.9999$ FID1939291959 $0.8-3.2$ $y = 1.32e + 06x - 5727$ $0.9999$ EthylbenzeneBID6132025140462 $0.4-2.3$ $y = 4.60e + 06x + 2456$ $0.9999$ FID1464803691216 $0.9-1.6$ $y = 1.34e + 06x - 5225$ $0.9999$ ButylbenzeneBID6132025140465 $0.3-3.3$ $y = 5.86e + 06x + 6786$ $0.9999$ FID2488203481139<	Decane							•	
FID78258261860 $0.3-17$ $y = 1.69e + 06 x - 12695$ $0.9994$ CyclohexaneBID637210379260 $0.4-3.7$ $y = 3.99e + 06 x - 16656$ $0.9993$ FID150496228753 $0.9-2.3$ $y = 6.48950 x - 3228$ $0.9993$ CyclodecaneBID6572168129426 $0.3-2.3$ $y = 6.04e + 06 x - 10544$ $0.9999$ CycloddecaneBID47515674191384 $0.2-7.7$ $y = 971333 x - 7485$ $0.9993$ CycloddecaneBID130429105346 $0.2-34$ $y = 6.33e + 06 x - 2164$ $0.9999$ BenzeneBID633208897320 $0.1-5.0$ $y = 4.06e + 06 x - 4824$ $0.9999$ FID25038261183603 $0.6-3.3$ $y = 1.26e + 06 x - 4874$ $0.9999$ TolueneBID4771574124408 $0.1-10$ $y = 4.40e + 06 x + 1766$ $0.9999$ FID119392291959 $0.8-3.2$ $y = 1.32e + 06 x - 5727$ $0.9999$ EthylbenzeneBID6132025140462 $0.4-2.3$ $y = 4.60e + 06 x - 5295$ $0.9999$ ButylbenzeneBID85728281139 $0.6-3.8$ $y = 1.69e + 06 x - 5377$ $0.9999$ PentylbenzeneBID7492472169559 $0.3-2.9$ $y = 5.72e + 06 x + 8451$ $0.9999$ PutylbenzeneBID7492472169<	~ .							2	
Cyclohexane   BID   637   2103   79   260   0.4–3.7   y = 3.99e + 06 x - 16656   0.9993     FID   150   496   228   753   0.9–2.3   y = 648950 x - 3228   0.9993     Cyclodecane   BID   657   2168   129   426   0.3–2.3   y = 6.04e + 06 x - 10544   0.9999     FID   475   1567   419   1384   0.2–7.7   y = 971333 x - 7485   0.9993     Cyclododecane   BID   130   429   105   346   0.2–34   y = 6.38e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2.3–6.2   y = 1.45e + 06 x - 48345   0.9999     Benzene   BID   633   2088   97   320   0.1–5.0   y = 4.06e + 06 x - 4874   0.9999     Toluene   BID   477   1574   124   408   0.1–10   y = 4.40e + 06 x + 1766   0.9999     FID   119   392   291   959   0.8–3.2   y = 1.32e + 06 x - 5225   0	Dodecane							2	
FID   150   496   228   753   0.9-2.3   y = 648950 x - 3228   0.9993     Cyclodecane   BID   657   2168   129   426   0.3-2.3   y = 6.04e + 06 x - 10544   0.9999     FID   475   1567   419   1384   0.2-7.7   y = 971333 x - 7485   0.9993     Cyclododecane   BID   130   429   105   346   0.2-34   y = 6.38e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2.3-6.2   y = 1.45e + 06 x - 48365   0.9999     Benzene   BID   633   2088   97   320   0.1-5.0   y = 4.06e + 06 x - 492   0.9999     FID   2503   8261   183   603   0.6-3.3   y = 1.26e + 06 x - 4874   0.9999     Toluene   BID   477   1574   124   408   0.1-10   y = 4.40e + 06 x + 1766   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4-2.3   y = 4.60e + 06 x + 5245 <th< th=""><th>Cucloboxana</th><th></th><th></th><th></th><th></th><th></th><th></th><th></th><th></th></th<>	Cucloboxana								
Cyclodecane   BID   657   2168   129   426   0.3-2.3   y = 6.04e + 06 x - 10544   0.9999     FID   475   1567   419   1384   0.2-7.7   y = 971333 x - 7485   0.9993     Cyclododecane   BID   130   429   105   346   0.2-34   y = 6.33e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2.3-6.2   y = 1.45e + 06 x - 4835   0.9994     Benzene   BID   633   2088   97   320   0.1-5.0   y = 4.06e + 06 x - 4874   0.9999     Toluene   BID   477   1574   124   408   0.1-10   y = 4.40e + 06 x + 1766   0.9999     Toluene   BID   477   1574   124   408   0.1-10   y = 4.40e + 06 x - 5227   0.9999     FID   119   392   291   959   0.8-3.2   y = 1.32e + 06 x - 5275   0.9999     BID   613   2025   140   462   0.4-2.3   y = 4.60e + 06 x + 45766   0	Cyclollexalle							•	
FID   475   1567   419   1384   0.2-7.7   y = 971333 x - 7485   0.9993     Cyclododecane   BID   130   429   105   346   0.2-34   y = 6.33e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2.3-6.2   y = 1.45e + 06 x - 8385   0.9994     Benzene   BID   633   2088   97   320   0.1-5.0   y = 4.06e + 06 x - 4874   0.9999     FiD   2503   8261   183   603   0.6-3.3   y = 1.26e + 06 x - 4874   0.9999     Toluene   BID   477   1574   124   408   0.1-10   y = 4.40e + 06 x + 1766   0.9999     FID   119   392   291   959   0.8-3.2   y = 1.32e + 06 x - 5727   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4-2.3   y = 4.60e + 06 x + 2456   0.9999     Butylbenzene   BID   857   2828   141   465   0.3-3.3   y = 5.86e + 06 x + 6786	Cyclodecane								
Cyclododecane   BID   130   429   105   346   0.2–34   y = 6.33e + 06 x - 2164   0.9999     FID   3634   11994   295   972   2.3–6.2   y = 1.45e + 06 x - 8385   0.9994     Benzene   BID   633   2088   97   320   0.1–5.0   y = 4.06e + 06 x - 492   0.9999     Toluene   BID   477   1574   124   408   0.1–10   y = 4.40e + 06 x + 1766   0.9999     Toluene   BID   477   1574   124   408   0.1–10   y = 4.40e + 06 x + 1766   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4–2.3   y = 4.60e + 06 x + 2456   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4–2.3   y = 4.60e + 06 x + 5295   0.9999     Butylbenzene   BID   857   2828   141   465   0.3–3.3   y = 5.86e + 06 x + 6786   0.9999     FID   219   723   345   1139   0.6–3	cyclouecalle								
FID3634119942959722.3-6.2 $y = 1.45e + 06 x - 8385$ 0.9994BenzeneBID6332088973200.1-5.0 $y = 4.06e + 06 x - 492$ 0.9999FID250382611836030.6-3.3 $y = 1.26e + 06 x - 4874$ 0.9999TolueneBID47715741244080.1-10 $y = 4.40e + 06 x + 1766$ 0.9999EthylbenzeneBID61320251404620.4-2.3 $y = 1.32e + 06 x - 5727$ 0.9999ButylbenzeneBID61320251404620.4-2.3 $y = 4.60e + 06 x + 2456$ 0.9999ButylbenzeneBID85728281414650.3-3.3 $y = 5.86e + 06 x + 6786$ 0.9999PentylbenzeneBID74924721695590.3-2.9 $y = 5.72e + 06 x - 5845$ 0.9999PentylbenzeneBID74924721695590.3-2.9 $y = 5.72e + 06 x + 8451$ 0.9999FID24882034811501.2-3.6 $y = 1.67e + 06 x - 5377$ 0.9999Ethyl acetateBIDn.d.n.d.42714080.5-1.2 $y = 7.58e + 06 x + 145864$ 0.998Butyl acetateBIDn.d.n.d.2113900.5-11 $y = 1.30e + 06 x - 7674$ 0.99992-Methoxy-1-propylBIDn.d.n.d.792620.5-3.5 $y = 6.61e + 06 x + 39141$ 0.9991	Cvclododecane								
Benzene   BID   633   2088   97   320   0.1-5.0   y = 4.06e + 06 x - 492   0.9999     FID   2503   8261   183   603   0.6-3.3   y = 1.26e + 06 x - 4874   0.9999     Toluene   BID   477   1574   124   408   0.1-10   y = 4.40e + 06 x + 1766   0.9999     FID   119   392   291   959   0.8-3.2   y = 1.32e + 06 x - 5727   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4-2.3   y = 4.60e + 06 x + 2456   0.9999     Butylbenzene   BID   857   2828   141   465   0.3-3.3   y = 5.86e + 06 x - 5295   0.9999     FID   146   480   369   1216   0.9-1.6   y = 1.32e + 06 x - 5845   0.9999     Butylbenzene   BID   857   2828   141   465   0.3-3.3   y = 5.86e + 06 x + 6786   0.9999     Pentylbenzene   BID   749   2472   169   559   0.3-2.9   y =	-,								
TolueneBID4771574124408 $0.1-10$ $y = 4.40e + 06 x + 1766$ $0.9999$ FID119392291959 $0.8-3.2$ $y = 1.32e + 06 x - 5727$ $0.9999$ EthylbenzeneBID6132025140462 $0.4-2.3$ $y = 4.60e + 06 x + 2456$ $0.9999$ ButylbenzeneBID1464803691216 $0.9-1.6$ $y = 1.34e + 06 x - 5295$ $0.9999$ ButylbenzeneBID8572828141465 $0.3-3.3$ $y = 5.86e + 06 x + 6786$ $0.9999$ PentylbenzeneBID7492472169559 $0.3-2.9$ $y = 5.72e + 06 x + 8451$ $0.9999$ FID2488203481150 $1.2-3.6$ $y = 1.67e + 06 x - 5377$ $0.9999$ Ethyl acetateBIDn.d.n.d.4471480 $0.5-1.2$ $y = 740187 x - 2275$ $0.9999$ Butyl acetateBIDn.d.n.d.102338 $0.6-6.5$ $y = 7.58e + 06 x + 37757$ $0.9991$ FIDn.d.n.d.79262 $0.5-3.5$ $y = 6.61e + 06 x + 39141$ $0.9991$	Benzene	BID		2088	97	320	0.1-5.0	y = 4.06e + 06 x - 492	
FID   119   392   291   959   0.8–3.2   y = 1.32e + 06 x - 5727   0.9999     Ethylbenzene   BID   613   2025   140   462   0.4–2.3   y = 4.60e + 06 x + 2456   0.9999     Butylbenzene   BID   146   480   369   1216   0.9–1.6   y = 1.34e + 06 x - 5295   0.9999     Butylbenzene   BID   857   2828   141   465   0.3–3.3   y = 5.86e + 06 x + 6786   0.9999     Pentylbenzene   BID   219   723   345   1139   0.6–3.8   y = 1.69e + 06 x - 5845   0.9999     Pentylbenzene   BID   749   2472   169   559   0.3–2.9   y = 5.72e + 06 x + 8451   0.9999     FID   248   820   348   1150   1.2–3.6   y = 1.67e + 06 x - 5377   0.9999     Ethyl acetate   BID   n.d.   n.d.   447   145   0.8–3.6   y = 6.31e + 06 x + 145864   0.998     FID   n.d.   n.d.   102   338   0.6–6.5		FID	2503	8261	183	603	0.6-3.3	y = 1.26e + 06 x - 4874	0.9999
Ethylbenzene   BID   613   2025   140   462   0.4–2.3   y = 4.60e + 06 x + 2456   0.9999     Butylbenzene   FID   146   480   369   1216   0.9–1.6   y = 1.34e + 06 x - 5295   0.9999     Butylbenzene   BID   857   2828   141   465   0.3–3.3   y = 5.86e + 06 x + 6786   0.9999     Pentylbenzene   BID   749   2472   169   559   0.3–2.9   y = 5.72e + 06 x + 8451   0.9999     FlD   248   820   348   1150   1.2–3.6   y = 1.67e + 06 x - 5377   0.9999     Ethyl acetate   BID   n.d.   n.d.   44   145   0.8–3.6   y = 6.31e + 06 x + 145864   0.998     Butyl acetate   BID   n.d.   n.d.   427   1408   0.5–1.2   y = 740187 x - 2275   0.9999     Butyl acetate   BID   n.d.   n.d.   102   338   0.6–6.5   y = 7.58e + 06 x + 37757   0.9991     Butyl acetate   BID   n.d.   n.d.	Toluene	BID	477	1574	124	408	0.1–10	y = 4.40e + 06 x + 1766	0.9999
FID   146   480   369   1216   0.9-1.6   y = 1.34e + 06 x - 5295   0.9999     Butylbenzene   BID   857   2828   141   465   0.3-3.3   y = 5.86e + 06 x + 6786   0.9999     FID   219   723   345   1139   0.6-3.8   y = 1.69e + 06 x - 5845   0.9999     Pentylbenzene   BID   749   2472   169   559   0.3-2.9   y = 5.72e + 06 x + 8451   0.9999     FID   248   820   348   1150   1.2-3.6   y = 1.67e + 06 x - 5377   0.9999     Ethyl acetate   BID   n.d.   n.d.   44   145   0.8-3.6   y = 6.31e + 06 x + 145864   0.998     Butyl acetate   BID   n.d.   n.d.   427   1408   0.5-1.2   y = 740187 x - 2275   0.9999     Butyl acetate   BID   n.d.   n.d.   102   338   0.6-6.5   y = 7.58e + 06 x + 37757   0.9991     Butyl acetate   BID   n.d.   n.d.   21   1390   0.5-11 <th></th> <th></th> <th>119</th> <th>392</th> <th>291</th> <th>959</th> <th></th> <th>y = 1.32e + 06 x - 5727</th> <th>0.9999</th>			119	392	291	959		y = 1.32e + 06 x - 5727	0.9999
Butylbenzene   BID   857   2828   141   465   0.3-3.3   y = 5.86e + 06 x + 6786   0.9999     FID   219   723   345   1139   0.6-3.8   y = 1.69e + 06 x - 5845   0.9999     Pentylbenzene   BID   749   2472   169   559   0.3-2.9   y = 5.72e + 06 x + 8451   0.9999     FID   248   820   348   1150   1.2-3.6   y = 1.67e + 06 x - 5377   0.9999     Ethyl acetate   BID   n.d.   n.d.   44   145   0.8-3.6   y = 6.31e + 06 x + 145864   0.998     FID   n.d.   n.d.   427   1408   0.5-1.2   y = 7.58e + 06 x + 37757   0.9999     Butyl acetate   BID   n.d.   n.d.   102   338   0.6-6.5   y = 7.58e + 06 x + 37757   0.9991     FID   n.d.   n.d.   421   1390   0.5-11   y = 1.30e + 06 x - 7674   0.9999     2-Methoxy-1-propyl   BID   n.d.   n.d.   79   262   0.5-3.5   y = 6.	Ethylbenzene								
FID2197233451139 $0.6-3.8$ $y = 1.69e + 06x - 5845$ $0.9999$ PentylbenzeneBID7492472169559 $0.3-2.9$ $y = 5.72e + 06x + 8451$ $0.9999$ FID2488203481150 $1.2-3.6$ $y = 1.67e + 06x - 5377$ $0.9999$ Ethyl acetateBIDn.d.n.d.44145 $0.8-3.6$ $y = 6.31e + 06x + 145864$ $0.998$ Butyl acetateBIDn.d.n.d.4271408 $0.5-1.2$ $y = 740187x - 2275$ $0.9999$ Butyl acetateBIDn.d.n.d.102338 $0.6-6.5$ $y = 7.58e + 06x + 37757$ $0.9991$ FIDn.d.n.d.4211390 $0.5-11$ $y = 1.30e + 06x - 7674$ $0.9999$ 2-Methoxy-1-propylBIDn.d.n.d. $79$ $262$ $0.5-3.5$ $y = 6.61e + 06x + 39141$ $0.9991$									
PentylbenzeneBID FID749 2482472 820169 348559 559 1.2-3.6 $y = 5.72e + 06 x + 8451$ $y = 1.67e + 06 x - 5377$ $y = 1.67e + 06 x - 5377$ 0.9999 0.9999Ethyl acetateBID FIDn.d.n.d.44 44145 1450.8-3.6 0.5-1.2 $y = 6.31e + 06 x + 145864$ $y = 740187 x - 2275$ 0.9999 0.9999Butyl acetateBID FIDn.d.n.d.102 421338 13900.6-6.5 0.5-11 $y = 740187 x - 2275$ $y = 7.58e + 06 x + 37757$ 0.9991 0.9991Pentyl acetateBID FIDn.d.n.d.102 421338 13900.5-11 0.5-11 $y = 1.30e + 06 x - 7674$ $y = 6.61e + 06 x + 391410.9991$	Butylbenzene							•	
FID thyl acetate248 BID FID820 n.d.348 n.d.1150 441.2-3.6 0.8-3.6 $y = 1.67e + 06x - 5377$ $y = 6.31e + 06x + 145864$ $y = 6.31e + 06x + 145864$ 0.9999 0.9999Butyl acetateBID FIDn.d.n.d.427 4271408 14080.5-1.2 0.5-1.2 $y = 740187x - 2275$ $y = 7.58e + 06x + 37757$ 0.9999 0.9999Butyl acetateBID FIDn.d.n.d.102 421338 13900.6-6.5 0.5-11 $y = 7.58e + 06x + 37757$ 0.9991 0.99912-Methoxy-1-propylBIDn.d.n.d.79 262262 0.5-3.5 $y = 6.61e + 06x + 39141$ 0.9991	Dentulhan								
Ethyl acetate   BID   n.d.   n.d.   44   145   0.8–3.6   y = 6.31e + 06 x + 145864   0.998     FID   n.d.   n.d.   427   1408   0.5–1.2   y = 740187 x - 2275   0.9999     Butyl acetate   BID   n.d.   n.d.   102   338   0.6–6.5   y = 7.58e + 06 x + 37757   0.9991     FID   n.d.   n.d.   421   1390   0.5–11   y = 1.30e + 06 x - 7674   0.9999     2-Methoxy-1-propyl   BID   n.d.   n.d.   79   262   0.5–3.5   y = 6.61e + 06 x + 39141   0.9991	rentyidenzene								
FID   n.d.   n.d.   427   1408   0.5-1.2   y = 740187 x - 2275   0.9999     Butyl acetate   BID   n.d.   n.d.   102   338   0.6-6.5   y = 7.58e + 06 x + 37757   0.9991     FID   n.d.   n.d.   121   1390   0.5-11   y = 1.30e + 06 x - 7674   0.9999     2-Methoxy-1-propyl   BID   n.d.   n.d.   79   262   0.5-3.5   y = 6.61e + 06 x + 39141   0.9991	Ethyl acotato								
Butyl acetate   BID   n.d.   n.d.   102   338   0.6–6.5   y = 7.58e + 06 x + 37757   0.9991     FID   n.d.   n.d.   421   1390   0.5–11   y = 1.30e + 06 x - 7674   0.9999     2-Methoxy-1-propyl   BID   n.d.   n.d.   79   262   0.5–3.5   y = 6.61e + 06 x + 39141   0.9991	Luiyi aceidle							-	
FIDn.d.n.d.4211390 $0.5-11$ $y = 1.30e + 06x - 7674$ 0.9999 <b>2-Methoxy-1-propyl</b> BIDn.d.n.d.79262 $0.5-3.5$ $y = 6.61e + 06x + 39141$ 0.9991	Butyl acetate								
<b>2-Methoxy-1-propyl</b> BID n.d. n.d. 79 262 0.5–3.5 $y = 6.61e + 06x + 39141$ 0.9991	baly acciate							2	
	2-Methoxy-1-propvl								
	acetate							,	

(continued)

Table 2. Continued.	Tab	le 2.	Continued.
---------------------	-----	-------	------------

		3ơ/S N	/lethod	S/N Method				
Analyte	Detector	LOD (pg/s)	LOQ (pg/s)	LOD (pg/s)	LOQ (pg/s)	RSD (%) range	Regression equation	Coefficient of determination R <sup>2</sup>
	FID	n.d.	n.d.	777	2565	0.5-1.1	$y = 908912 \ x - 5950$	0.9999
2-Butoxyethyl acetate	BID	n.d.	n.d.	150	494	0.4–3.4	y = 6.60e + 06 x + 49378	0.999
	FID	n.d.	n.d.	839	2767	0.1-8.7	y = 1.14e + 06 x - 6964	0.9999

LOD: Limit of detection

LOQ: Limit of quantification

RSD: Relative standard deviation

BID: Barrier discharge ionization detector

FID: Flame ionization detector


>0.9999 for anilines, halogenated compounds, phenols, alcohols, cyclic, and aromatic compounds with the BID and for four compounds also with the FID.

The relative standard deviation (RSD) was determined and the results were found acceptable for both detectors with the majority of the RSD values being lower than 5%, except for the measurements at the lowest concentrations injected  $(1 \text{ ng}/\mu\text{L}, 100 \text{ pg}/\mu\text{L})$ , where the relative precision increased to above 10% (Table 2). Particularly, the halogenated compounds and aldehydes had large standard deviations for repeated measurements at low concentrations. A large part of the encountered variability is due to the low concentration levels chosen for these measurements, considering that already for statistical reasons, the RSD approaches the value of 33% when the concentration is close to the limit of quantification (LOQ), and may rise even to 50% at the lower limit of the guantitative working range. The RSD values for most analytes were lower for the BID. For the standards measured at two concentrations, the RSD was also investigated after ten days and it changed less than 20% for the 1 ng/ $\mu$ L concentration.


Figure 6 compares the sensitivity of the BID to that of the FID for the standards at five concentration levels. In general, the sensitivity of the BID is higher than the FID sensitivity. The sensitivity was determined as the slope of the calibration curve, and the results vary depending on the analyte. The BID provides for all analytes and analyte classes a higher response than the FID, although the difference in sensitivity depends on the particular class of compounds. Thus, the sensitivity of the BID detector was found to be higher than the sensitivity of the FID by a factor of 3.1 for anilines, 4.1 for aldehydes, 3.7 for halogens, 3.8 for phenols, 6.4 for alcohols, 4.1 for alkanes, 4.1 for cyclic compounds, 2.4 for aromatics, and 5.9 for esters. Generalizing this finding over all classes of compounds, the sensitivity of the BID was higher on average by a factor of 4.2 compared to the FID.

In the next step, response factors (RFs, determined as peak area/mass) were calculated for the individual compounds and investigated more closely to detect any systematic trends. Figure 7 reports the RFs for four analytes from different chemical classes. Although not completely homogeneous, no concentration dependence of the RFs can be observed, and the variations observed are more likely an effect of the chromatographic conditions, rather than the detection. The obtained RFs for the BID were generally higher, reflecting the higher sensitivity of the BID versus the FID.


Using the results from the measurement of standards at a higher concentration level (100 ng/ $\mu$ L), the difference in RFs for both detectors is depicted in Figure 8 in form of the



**Figure 6.** Normalized sensitivity of the BID in comparison to the FID. Normalized sensitivity: (compound sensitivity-lowest sensitivity from all 30 compounds)/lowest sensitivity from all 30 compounds. Compound numbers: (1) N,N-dimethylaniline, (2) N,N-diethylaniline, (3) isovalevaldehyde, (4) benzaldehyde, (5) 1-bromo-3-chloropropane, (6) 1-bromohexane, (7) 1-bromodecane, (8) phenol, (9) 2,4-dichlorophenol, (10) 4-chlorophenol, (11) 2-propanol, (12) 1-butanol, (13) isoamyl alcohol, (14) 1,4-butanediol, (15) 1-octanol, (16) heptane, (17) decane, (18) dodecane, (19) cyclohexane, (20) cyclodecane, (21) cyclododecane, (22) benzene, (23) toluene, (24) ethylbenzene, (25) butylbenzene, (26) pentylbenzene, (27) ethyl acetate, (28) butyl acetate, (29) 2-methoxy-1-propyl acetate, and (30) 2-butoxyethyl acetate.



**Figure 7.** Concentration dependence of the response factors (peak area/mass) for four compounds at five concentrations. The error bars represent  $\pm$  one standard deviation.



response factor ratio RF(BID)/RF(FID). In Figure 8(c), the RFs of different classes of chemical compounds are reported. It is evident that the BID RFs, although generally higher than those of the FID, differ significantly. The compound class that has the highest response ratio, and thus the highest relative sensitivity in the comparison of BID and FID, is the one of the alcohols, followed by the halogenated compounds. Somewhat surprisingly, the PAHs had the lowest and most constant RF ratio between BID and FID.

The above observations can be explained through the different ionization mechanisms in the BID and the FID. Since the main route of ionization in the BID is through photoionization, it can be expected that a constant number of charges (cations and electrons) are produced per mole of analyte molecules. In contrast to this, the FID generates its signal by completely atomizing a molecule and simultaneously ionizing the carbon fragments produced. The FID response is thus not concentration-dependent (mol/s), but mass-dependent (g/s).


Moreover, the FID response is also diminished compared to pure hydrocarbons for analytes that contain heteroelements. According to our chemical understanding, a carbon atom bound to one or more heteroatoms such as O, S, Cl, or Br is considered already as oxidized. It cannot be further oxidized, which reduces incrementally the response of the FID. This concept is already known for long under the term "effective carbon number" (ECN) concept (Scanlon and Willis 1985) and can be used to theoretically calculate the FID RF of an analyte as a function of the carbon number and the functional groups attached to the carbon backbone.

The BID/FID RF ratios show remarkable trends within the individual compound classes: For example, the RF ratio clearly increases for the chlorinated compounds investigated with the degree of substitution by chlorine while it decreases for some of the other compound classes with chain length.

In order to understand which phenomenon is responsible for such behavior, Figure 8 presents the original RFs of the BID (Figure 8(a) and the FID (Figure 8(b)). For the halogenated compounds—which in the current case are mono- to tetrasubstituted benzenes—this is explained by a strongly decreasing RF of the FID with increase in the chlorine substitution of the aromatic ring.

A contrary effect is seen for the aliphatic alcohols and the *n*-alkanes whose response decreases with increasing chain length. Since this can be caused both by a decrease of the response of the BID or by the increase of the RF of the FID, the RFs of the two instruments were evaluated and plotted separately. The results show that the main influence on this trend is the significant change in specific response for the BID. This change is due to the fact that the RFs were normalized to *mass* of analyte, *not* to *amount of substance*. Given that the BID leads to the formation of single-charged cations, the amount of ions generated

**Figure 8.** Diagram with the specific peak area (area per mass) values per compound for both detectors. Analytes: (1) cyclohexane, (2) cycloheptane, (3) cyclododecane, (4) cyclopentanone, (5) cyclohexanone, (6) cyclodecanone, (7) cyclododecanone, (8) chlorobenzene, (9) 1,2-dichlorobenzene, (10) 1,2,3-trichlorobenzene, (11) 1,2,3,4-tetrachlorobenzene, (12) chloroheptane, (13) chlorooctane, (14) chlorodecane, (15) chlorododecane, (16) 1-pentanol, (17) 1-hexanol, (18) 1-heptanol, (19) 1-octanol, (20) 1-nonanol, (21) 1-decanol, (22) n-hexane, (23) n-heptane, (24) *n*-octane, (25) *n*-decane, (26) n-dodecane, (27) n-tetradecane, (28) naphthalene, (29) acenaphthene, (30) fluorene, (31) phenanthrene, (32) fluoranthene, (33) pyrene, (34) chrysene, (35) benzene, (36) toluene, (37) ethylbenzene, (38) propylbenzene, (39) butylbenzene, and (40) pentylbenzene.



**Figure 9.** Response of the barrier discharge ionization detector (BID) and flame ionization detector (FID) for alkane standards C10-C17 at  $1 \text{ ng}/\mu \text{L}$ . The analytical conditions are the same as in Figure 5.


is proportional to the amount of substance of this analyte, but inversely proportional to the molecular mass of this compound. The longer the aliphatic chain of the analyte, the lower its specific RF. This appears to be the main reason for the observed decrease in RF ratios for the *n*-alcohols and *n*-alkanes with increasing molecular mass or chain length.

A further experiment was performed in order to verify that the BID responds to the amount of substance arriving at the detector per time unit, rather than mass of analyte per time unit. To this end, an alkane standard with the same mass concentration and carbon numbers from ten to seventeen was analyzed by both detectors. Figure 9 exhibits a decrease of the response with increasing carbon number in the BID in contrast to the FID response where this decrease is not evident. The reason for this difference is that the BID, in contrast to the FID, does not lead to the complete atomization of the analytes because of its low thermal energy and excitation. Thus, only one cation is formed per molecule in the BID with one electron being produced in the same process while a  $C_n$  compound can produce n carbocations and n electrons in the BID. The ionization yield is clearly lower than unity but can be expected to be similar for both excitation sources.

#### Limits of detection and quantification

The limit of detection (LOD) was determined by two different approaches. The first was based on the definition of the LOD as the concentration at which the signal-tonoise ratio (S/N) is equal to 3. Both signal and noise are expressed in peak height units. The second approach uses measurements of the standard deviation of the blank  $\sigma$  (n=5) to calculate the LOD as  $3\sigma/S$  (where S is slope or sensitivity) from the regression equation (Figure 10; Table 2). Also, in this case, the regression line was calculated for the peak height.

Given that the calculation of the LOD was performed in two different ways, there are evident differences between the results. According to the LOD values obtained with the S/N ratio method, the BID always exhibited better detection limits than the FID, especially for the halogenated, oxygenated compounds, and for some alkanes where the BID LOD was one order of magnitude lower than for the FID. The LOD values determined



**Figure 10.** Comparison of the limits of detection (LODs) for the two detectors: (a) calculation when the value of the signal-to-noise ratio (S/N) equals to 3. The compounds on the *x* axis are identified in Figure 6. (b) Use of the height of the blank measurement (n = 5) to calculate the limit of detection using the  $3\sigma/S$  relationship where S is the slope or sensitivity.

from the blank covered two orders of magnitude, ranging from the medium pg/s to the low ng/s range.

The LOD in the  $3\sigma/S$  results is influenced by the background in the blank samples. The background in the BID blanks was high as a consequence of residual impurities in the He carrier and plasma gas and hence the LOD values were higher. For the FID, the background was low in comparison to the BID and the LODs were lower in comparison to the values obtained with the other method. Generally, the FID values were lower with exceptions for compounds such as heptane and cyclododecane. Because of the high background in the blanks, the LOD values do not reflect the same increase in sensitivity of the BID compared to the FID as when evaluating the slopes of the calibration lines of the two detectors which is typically used to express the instrumental or method sensitivity.

#### Conclusions

The current study aimed to provide a better understanding and a more comprehensive and quantitative characterization of the response behavior of the BID and to compare its performance to that of the FID for gas chromatography.

The detection principle of the BID was identified as a photoionization process of the analyte molecules, caused by metastable He species created in the dielectric barrier discharge of the detector unit. As a result, the BID shows a concentration-dependent rather than mass-dependent response behavior compared to the FID. In general, the BID shows a sensitivity that is about four times higher than that of the FID, which is particularly the case for the compounds to which the FID has a medium to low response, notably oxygen- and halogen-containing compounds.

The BID response also requires optimization and is dependent on various factors, such as the detector temperature and detector make-up gas flow rate. Depending on the particular compound, detection limits for the BID in the two-digit picogram range (absolute) were achieved. We expect that limits of detection in the low one-digit pg (absolute) range can be reached. However, to this end, even higher requirements are imposed to the carrier gas supply—the use of He with purity of better than 99.9999% rather than better than 99.9999% as is often recommended for GC-FID is—and to the stability of the chromatographic columns (i.e., absence of bleed).

In general, the BID has similar robustness and ease of operation as the FID. Two practical advantages of the BID compared to the FID are the absence of a burner head which may clog with time, and that only He is required as detector gas which is also the GC carrier gas. The investment for the commercial BID instrument is moderately higher than that for a regular FID.

To summarize, the BID is a highly versatile and an interesting addition to the toolkit of gas chromatographic detectors. Due to its functioning principles, it is best used for low-molecular weight, heteroatom containing analytes where it easily outperforms the FID, while this advantage is gradually lost with increasing molecular weight of the analyte. Moreover, achieving highest performance with the BID imposes even more stringent requirements and control of the operating conditions and system state than for GC-FID. Where the above stated conditions are met—e.g., in the analysis of low molecular weight ozone precursors—the BID will find highly useful niches of application.

#### Acknowledgements

One of the authors (M.A.) would like to acknowledge an Erasmus fellowship that allowed her to carry out this work. Part of this study was carried out within the Austrian Research Promotion Agency (FFG) funded project "DianaBatt" (Project number 858298) which is also gratefully acknowledged.

#### ORCID

Maria Antoniadou (b) http://orcid.org/0000-0003-0416-5120 George A. Zachariadis (b) http://orcid.org/0000-0001-6191-6301 Erwin Rosenberg (b) http://orcid.org/0000-0002-3608-4571 2838 🛞 M. ANTONIADOU ET AL.

#### References

- Brandt, S., A. Schütz, F. D. Klute, J. Kratzer, and J. Franzke. 2016. Dielectric barrier discharges applied for optical spectrometry. *Spectrochimica Acta Part B* 123:6–32. doi: 10.1016/j.sab.2016.07.001.
- Franchina, F. A., M. Maimone, D. Sciarrone, G. Purcaro, P. Q. Tranchida, and L. Mondello. 2015. Evaluation of a novel helium ionization detector within the context of (low-)flow modulation comprehensive two-dimensional gas chromatography. *Journal of Chromatography* A 1402:102–9. doi: 10.1016/j.chroma.2015.05.013.
- Gras, R., J. Luong, M. Monagle, and B. Winniford. 2006. Gas chromatographic applications with the dielectric barrier discharge detector. *Journal of Chromatographic Science* 44(2):101–7. doi: 10.1093/chromsci/44.2.101.
- Guo, C., F. Tang, J. Chen, X. Wang, S. Zhang, and X. Zhang. 2015. Development of dielectricbarrier-discharge ionization. *Analytical and Bioanalytical Chemistry* 407(9):2345–64. doi: 10.1007/s00216-014-8281-y.
- Hagenhoff, S., and H. Hayen. 2018. LC/MS analysis of vitamin D metabolites by dielectric barrier discharge ionization and a comparison with electrospray ionization and atmospheric pressure chemical ionization. *Analytical and Bioanalytical Chemistry* 410(20):4905–11. doi: 10.1007/s00216-018-1137-0.
- Han, B., X. Jiang, X. Hou, and C. Zheng. 2014. Dielectric barrier discharge carbon atomic emission spectrometer: universal GC detector for volatile carbon-containing compounds. *Analytical Chemistry* 86(1):936–42. doi: 10.1021/ac403662w.
- Hu, J., W. Li, C. Zheng, and X. Hou. 2011. Dielectric barrier discharge in analytical spectrometry. *Applied Spectroscopy* 46(5):368–87. doi: 10.1080/05704928.2011.561511.
- Huba, A. K., M. F. Mirabelli, and R. Zenobi. 2018. High-throughput screening of PAHs and polar trace contaminants in water matrices by direct solid-phase microextraction coupled to a dielectric barrier discharge ionization source. *Analytica Chimica Acta* 1030:125–32. doi: 10.1016/j.aca.2018.05.050.
- Iguchi, S., K. Teramura, S. Hosokawa, and T. Tanaka. 2015. Effect of the chloride ion as a hole scavenger on the photocatalytic conversion of CO2 in an aqueous solution over Ni–Al layered double hydroxides. *Physical Chemistry Chemical Physics* 17(27):17995–8003. doi: 10.1039/C5CP02724A.
- Jiang, X., C. Li, Z. Long, and X. Hou. 2015. Selectively enhanced molecular emission spectra of benzene, toluene and xylene with nano-MnO<sub>2</sub> in atmospheric ambient temperature dielectric barrier discharge. Analytical Methods 7(2):400–4. doi: 10.1039/C4AY02194K.
- Jiang, X., Y. Chen, C. Zheng, and X. Hou. 2014. Electrothermal vaporization for universal liquid sample introduction to dielectric barrier discharge microplasma for portable atomic emission spectrometry. Analytical Chemistry 86(11):5220–4. doi: 10.1021/ac500637p.
- Jo, S.-H., and K.-H. Kim. 2017. The applicability of a large-volume injection (LVI) system for quantitative analysis of permanent gases  $O_2$  and  $N_2$  using a gas chromatograph/barrier discharge ionization detector. *Environmental Monitoring and Assessment* 189(7):317.
- Kirk, A. T., T. Last, and S. Zimmermann. 2017. A sensitive gas chromatography detector based on atmospheric pressure chemical ionization by a dielectric barrier discharge. *Journal of Chromatography A* 1483:120–6. doi: 10.1016/j.chroma.2016.12.071.
- Kogelschatz, U. 2003. Dielectric-barrier discharges: their history, discharge physics, and industrial applications. *Plasma Chemistry and Plasma Processing* 23(1):1–46.
- Kratzer, J., S. Musil, K. Marschner, M. Svoboda, T. Matoušek, Z. Mester, R. E. Sturgeon, and J. Dědina. 2018. Behavior of selenium hydride in heated quartz tube and dielectric barrier discharge atomizers. *Analytica Chimica Acta* 1028:11–21. doi: 10.1016/j.aca.2018.05.053.
- Li, C., X. Jiang, and X. Hou. 2015. Dielectric barrier discharge molecular emission spectrometer as gas chromatographic detector for amines. *Microchemical Journal* 119:108–13. doi: 10.1016/j.microc.2014.11.005.

- Li, W., C. Zheng, G. Fan, L. Tang, K. Xu, Y. Lv, and X. Hou. 2011. Dielectric barrier discharge molecular emission spectrometer as multichannel GC detector for halohydrocarbons. *Analytical Chemistry* 83(13):5050–5. doi: 10.1021/ac2007224.
- Liu, Z., Z. Zhu, Q. Wu, S. Hu, and H. Zheng. 2011. Dielectric barrier discharge-plasma induced vaporization and its application to the determination of mercury by atomic fluorescence spectrometry. *The Analyst* 136(21):4539–44. doi: 10.1039/c1an15332c.
- Martin, A. J. P., and R. L. M. Synge. 1941. A new form of chromatogram employing two liquid phases. *Biochemical Journal* 35(12):1358–68. doi: 10.1042/bj0351358.
- Pascale, R., M. Caivano, A. Buchicchio, I. M. Mancini, G. Bianco, and D. Caniani. 2017. Validation of an analytical method for simultaneous high-precision measurements of greenhouse gas emissions from wastewater treatment plants using a gas chromatography-barrier discharge detector system. *Journal of Chromatography A* 1480:62–9. doi: 10.1016/j.chroma.2016.11.024.
- Rosenberg, E. 2006. Trace organic analysis by gas chromatography with selective detectors. In *Encyclopedia of analytical chemistry*, eds. R. A. Meyers and M. P. Miller. Hoboken, NJ: Wiley.
- Scanlon, J. T., and D. E. Willis. 1985. Calculation of flame ionization detector relative response factors using the effective carbon number concept. *Journal of Chromatographic Science* 23(8): 333–40. doi: 10.1093/chromsci/23.8.333.
- Schütz, A., S. Brandt, S. Liedtke, D. Foest, U. Marggraf, and J. Franzke. 2015. Dielectric barrier discharge ionization of perfluorinated compounds. *Analytical Chemistry* 87(22):11415–9. doi: 10.1021/acs.analchem.5b03538.
- Shinada, K., S. Horiike, S. Uchiyama, R. Takechi, and T. Nishimoto. 2012. Development of new ionization detector for gas chromatography by applying dielectric barrier discharge. *Shimadzu Review* SR13\_001E:1–5.
- Siemens, W. 1857. Ueber die elektrostatische Induktion und die Verzögerung des Stroms in Flaschendrähten. Annalen der Physik und Chemie 178(9):66–122. doi: 10.1002/andp.18571780905.
- Straka, M., S. Burhenn, K. Marschner, S. Brandt, U. Marggraf, J. Dědina, J. Franzke, and J. Kratzer. 2018. Novel designs of dielectric barrier discharge hydride atomizers for atomic spectrometry. Spectrochimica Acta Part B 146:69–76. doi: 10.1016/j.sab.2018.05.006.
- Ueta, I., Y. Nakamura, K. Fujimura, S. Kawakubo, and Y. Saito. 2017. Determination of gaseous formic and acetic acids by a needle-type extraction device coupled to a gas chromatography-barrier discharge ionization detector. *Chromatographia* 80(1):151–6. doi: 10.1007/s10337-016-3201-2.
- Ueta, I., Y. Nakamura, S. Kawakubo, and Y. Saito. 2018. Determination of aqueous formic and acetic acids by purge-and-trap analysis with a needle-type extraction device and gas chromatography barrier discharge ionization detector. *Analytical Sciences* 34(2):201–5. doi: 10.2116/analsci.34.201.
- Weatherly, C. A., R. M. Woods, and D. W. Armstrong. 2014. Rapid analysis of ethanol and water in commercial products using ionic liquid capillary gas chromatography with thermal conductivity detection and/or barrier discharge ionization detection. *Journal of Agricultural and Food Chemistry* 62(8):1832–8. doi: 10.1021/jf4050167.
- Zhu, H., M. Zhou, J. Lee, R. Nidetz, K. Kurabayashi, and X. Fan. 2016. Low-power miniaturized helium dielectric barrier discharge photoionization detectors for highly sensitive vapor detection. *Analytical Chemistry* 88(17):8780–6. doi: 10.1021/acs.analchem.6b02180.