

Development of an indirect enzyme-linked immunosorbent assay for the detection of osmotolerant yeast *Zygosaccharomyces rouxii* in different food

Hong Guo, Ya-Hong Yuan, Chen Niu, Yue Qiu, JianPing Wei & Tian-Li Yue

To cite this article: Hong Guo, Ya-Hong Yuan, Chen Niu, Yue Qiu, JianPing Wei & Tian-Li Yue (2018) Development of an indirect enzyme-linked immunosorbent assay for the detection of osmotolerant yeast *Zygosaccharomyces rouxii* in different food, Food and Agricultural Immunology, 29:1, 976-988, DOI: [10.1080/09540105.2018.1491532](https://doi.org/10.1080/09540105.2018.1491532)

To link to this article: <https://doi.org/10.1080/09540105.2018.1491532>



© 2018 The Author(s). Published by Informa UK Limited, trading as Taylor & Francis Group



Published online: 08 Oct 2018.



Submit your article to this journal [↗](#)



Article views: 502



View related articles [↗](#)



View Crossmark data [↗](#)

Development of an indirect enzyme-linked immunosorbent assay for the detection of osmotolerant yeast *Zygosaccharomyces rouxii* in different food

Hong Guo^{a,b,c}, Ya-Hong Yuan^{a,b,c}, Chen Niu^{a,b,c}, Yue Qiu^{a,b,c}, JianPing Wei^{a,b,c} and Tian-Li Yue^{a,b,c}

^aCollege of Food Science and Engineering, Northwest A & F University, Yangling, People's Republic of China; ^bLaboratory of Quality & Safety Risk Assessment for Agro-products (Yangling), Ministry of Agriculture, Yangling, People's Republic of China; ^cNational Engineering Research Center of Agriculture Integration Test (Yangling), Yangling, People's Republic of China

ABSTRACT

A polyclonal antibody-based indirect ELISA for the sensitive and specific detection of osmotolerant *Zygosaccharomyces rouxii* was established. The reactivity of the polyclonal antibodies (pAbs) with *Z. rouxii* was confirmed by immunofluorescence. The ELISA detection limit was 5×10^4 CFU/mL. For an evaluation of the specificity of pAbs, reactivity to 10 fungi species and 2 bacteria species commonly found in food environments was analysed by indirect ELISA. ELISA reacts with all tested *Z. rouxii*, whereas no cross-reactions or weakly reactions with 9 fungi and 2 bacteria were observed. The immunoassay staining confirmed the ELISA. Further on, the assay applicability was tested with concentrated apple juice, wine, and honey. Food samples were inoculated with 100 CFU/mL of *Z. rouxii*. After enrichment for 24 h, the established ELISA can successfully detect contaminated wine and concentrated apple juice samples, while only 30 h was required for honey samples.

ARTICLE HISTORY



Received 24 March 2018
Accepted 23 May 2018


KEYWORDS

Zygosaccharomyces rouxii;
polyclonal antibodies;
indirect ELISA; immuno-
detection; food quality

1. Introduction

Zygosaccharomyces rouxii is known for its ability to survive in the face of stress caused by high concentrations of ionic (salt) and non-ionic (sugar and polyols) (Dakal, Solieri, & Giudici, 2014). Combined with its ability to grow at low pH conditions and preservatives, it can carry out fermentation thus increases the risk of juice industries or cause high sugary and salty foods spoilage. For instance, in the traditional fermented foods, *Z. rouxii* produce aromatic flavours and alcohol (Watanabe, Uehara, & Mogi, 2013). In addition, the flor

CONTACT Tian-Li Yue  yueltl305@nwsuaf.edu.cn, yt6503@163.com  College of Food Science and Engineering, Northwest A & F University, Yangling, Shaanxi 712100, People's Republic of China, Laboratory of Quality & Safety Risk Assessment for Agro-products (Yangling), Ministry of Agriculture, Yangling, Shaanxi 712100, People's Republic of China, National Engineering Research Center of Agriculture Integration Test (Yangling), Yangling, Shaanxi 712100, People's Republic of China

 Supplemental data for this article can be accessed at <http://dx.doi.org/10.1080/09540105.2018.1491532>.

© 2018 The Author(s). Published by Informa UK Limited, trading as Taylor & Francis Group

This is an Open Access article distributed under the terms of the Creative Commons Attribution-NonCommercial-NoDerivatives License (<http://creativecommons.org/licenses/by-nc-nd/4.0/>), which permits non-commercial re-use, distribution, and reproduction in any medium, provided the original work is properly cited, and is not altered, transformed, or built upon in any way.

formation spoilage caused by *Z. rouxii* in the surface of soy sauce (Watanabe et al., 2013) or wine (Zuehlke, Petrova, & Edwards, 2013) have been reported. More recently, the spoilage ability of *Z. rouxii* has been also investigated by our laboratory (Wang, Hu, Long, Niu, et al., 2015). The *Z. rouxii* but not the other 22 species of osmotolerant isolates could survive in concentrated apple juice. Although the spoilage of food industries caused by *Z. rouxii* can be noticed by swelling of packaging or flor formation, the products would have already been seriously spoiled at that point time, resulting in serious economic losses (Martorell, Stratford, Steels, Fernandez-Espinar, & Querol, 2007).

Therefore, it is essential for food manufacturers to implement routine inspection to detect potential *Z. rouxii* contamination in raw material, finished products and processing environment. Traditional microbiological technique used for quality control tests is cost-effective but laborious and time consuming, especially for *Z. rouxii* because it often takes a long time (48–72 h) to grow under high sugar conditions (Wang, Hu, Long, Guo, et al., 2015). And the culturing method has an obvious drawback that the actual number of microbial communities is usually underestimated (Gobbi et al., 2010). Researchers have made great efforts in the rapid identification of *Z. rouxii* and most reports rely on genetic-based methods. Pearson and McKee (1992) used two primers from *Z. rouxii* plasmid in a multiplex polymerase chain reaction (PCR) to identify this yeast. Identification of *Zygosaccharomyces* species including *Z. rouxii* species is normally done by sequencing the IST or D1/D2 rDNA regions (Harrison, Muir, Stratford, & Wheals, 2011; Hulin & Wheals, 2014). However, these genetic-based assays require skilled personnel and may not be easily adapted to routine testing. Other available detection methods, including gas chromatography–mass spectrometry (GC-MS) (Zhang, Yue, & Yuan, 2013), high-performance liquid chromatography (HPLC) (Arroyomanzanas, Huertaspérez, Gámizgracia, & Garcíacampaña, 2015), and electronic nose (Gobbi et al., 2015; Wang et al., 2016) are expensive and/or require well-trained laboratory technicians. Therefore, a more affordable, rapid, and sensitive method for detection of *Z. rouxii* is required.

Enzyme-linked immunosorbent assay (ELISA) have been used to detect food-borne microorganisms such as *Vibrio parahaemolyticus*, *Alicyclobacillus spp.*, *Escherichia coli* (Kumar et al., 2011; Li, Huang, Xia, & Liu, 2014; Li, Xia, & Yu, 2013; Wang, He, & Shi, 2007). Due to its sensitivity, specificity, rapidity, and simplicity, it may become a promising alternative of plating detection method. More recently, ELISA was applied to detect spores of *A. acidoterrestris* in complex food matrix (Mast, Dietrich, Didier, & Martlbauer, 2016). These ELISAs based on polyclonal, monoclonal antibodies or even antibody fragments can be well applied to the detection of bacteria in complex food matrices.

However, fungi share conserved glycoprotein immunogens, which means that the antibodies were often unable to differentiate the target fungus from non-target fungus (Thornton & Wills, 2015). The immuno-detection of fungi usually presents the cross-reaction problem between pAbs and targets (Schmechel, Green, Blachere, Janotka, & Beezhold, 2008; Vojdani, 2004). Hoffmann et al. (2016) and Wang et al. (2013) have suggested that structures of antigens or analytes would influence the selectivity of the antibodies, resulting cross-reactions. Maybe for this reason, there are very few reports focused on the pAbs-based immunoassays for food fungi testing. Here in, we used *Z. rouxii* B-WHX-12-54 grown under high sugar (600 g/L) stress as an immunogens to obtain rabbit pAbs. The reactivity of pAbs to 9 fungi and 2 bacteria commonly found in food environments was analysed by Indirect ELISA. We found that pAbs only obviously

recognise *Z. rouxii*, suggesting polyclonal antibody-based ELISA can be used in the detection of *Z. rouxii* in food. Since *Z. rouxii* leads to spoilage of sugar containing or even high sugar foods (apple juice concentrated, wine and honey) (Dakal et al., 2014; Loureiro & Malfeitoferreira, 2003), we established the indirect ELISA to detect *Z. rouxii* in enriched samples of these foods.

2. Materials and methods

2.1. Chemicals and reagents

Bovine serum albumin (BSA), Goat anti-rabbit IgG (H + L)-horseradish peroxidase (Goat anti-rabbit IgG (H + L)-HRP), and 3,3',5,5'-tetramethyl benzidine (TMB) were purchased from Sigma (St. Louis, MO, USA). The BCA protein assay kit and ECL western blotting detection kit were purchased from Pierce (Rockford, IL., USA). Tween-20, H₂SO₄ and commonly used reagents were purchased from local supplies at high purity. The coating buffer was Phosphate-buffered saline (PBS, pH 7.4). Phosphate-buffered saline (PBS, pH 7.4) with 0.05% (v/v) Tween-20 (PBST) was used for washing buffer and 2 mol/L H₂SO₄ was used as the halting solution. Concentrated apple juice (70° Brix) were purchased from a local juice factory (Shaanxi Haisheng Fresh Fruit Juice Co., Ltd.). Milk powder, Wine and honey were purchased from local supermarket.

2.2. Microorganisms and culture

All of the yeast strains and bacterial strains used are listed in Table 1. Fourteen wild-type isolates of different *Z. rouxii*, *Saccharomyces cerevisiae*, *Pichia kluyveri*, *P. kudriavzevii*, *Kluyveromyces marxianus*, *Candida tropicalis*, *C. orthopsilosis*, *Hanseniaspora uvarum*, *Debaryomyces hansenii* as well as *Wickerhamomyces anomalus* were obtained from Microbiology Laboratory of the College of Food Science and Engineering, Northwest A & F Univ. (NWSUAF) (Wang et al., 2012; Wang, Hu, Long, Niu, et al., 2015). Three strains of different *Z. rouxii* as well as two strains of *S. cerevisiae* and *E. coli* were purchased from the American Type Culture Collection (ATCC). From the German Collection of Microorganisms and Cell Cultures (DSMZ), single strain of *Alicyclobacillus acidoterrestris*

Table 1. Strains used in this study.

Species	Strains/isolates	Species	Strains/isolates
<i>Zygosaccharomyces rouxii</i>	ATCC ^a 2623	<i>Kluyveromyces marxianus</i>	B-WHX-12-21 ^b
<i>Zygosaccharomyces rouxii</i>	ATCC ^a 8383	<i>Pichia kudriavzevii</i>	B-WHX-12-36 ^b
<i>Zygosaccharomyces rouxii</i>	B-WHX-12-54 ^b	<i>C. tropicalis</i>	B-WHX-12-52 ^b
<i>Zygosaccharomyces rouxii</i>	B-WHX-12-53 ^b	<i>Hanseniaspora uvarum</i>	B-WHX-12-63 ^b
<i>Zygosaccharomyces rouxii</i>	B-WHX-12-55 ^b	<i>Hanseniaspora uvarum</i>	B-WHX-12-14 ^b
<i>Zygosaccharomyces rouxii</i>	ATCC ^a 66351	<i>Hanseniaspora uvarum</i>	B-WHX-12-61 ^b
<i>Saccharomyces cerevisiae</i>	1917 ^b	<i>Debaryomyces hansenii</i>	ACCC ^d 20310
<i>Saccharomyces cerevisiae</i>	B-WHX-12-05 ^b ATCC ^a 38531	<i>Wickerhamomyces anomalus</i>	B-WHX-12-11 ^b
<i>Saccharomyces cerevisiae</i>	B-WHX-12-16 ^b	<i>Escherichia coli</i>	ATCC ^a 25922
<i>Pichia kluyveri</i>	B-WHX-12-30	<i>Alicyclobacillus acidoterrestris</i>	DSM ^c 3924
<i>Candida. orthopsilosis</i>			

^aATCC: American Type Culture Collection.

^bNWSUAF: College of Food Science and Engineering, Northwest A&F Univ.

^cDSM: German Collection of Microorganisms and Cell Cultures.

^dACCC: Agricultural Culture Collection of China.

was acquired. Single strain of *D. hansenii* was purchased from Agricultural Culture Collection of China (ACCC). In brief, all yeast strains were inoculated in Yeast extract peptone dextrose broth (YPD, 1% yeast extract, 2% peptone and 2% glucose) liquid medium and cultured at 28 °C for 48 h. *A. acidoterrestris* was cultured in *A. acidoterrestris* medium (AAM, 4.0 g yeast extract, 2.0 g glucose, 0.4 g ammonia sulphate, 0.5 g calcium chloride, 1 g magnesium sulphate, 1.2 g monopotassium phosphate and 0.5 g manganese sulphate per litre of deionised water). *E. coli* was cultured in luria broth (LB, 1% tryptone, 0.5% yeast extract and 1% NaCl).

2.3. Polyclonal antibodies production

Z. rouxii B-WHX-12-54 was cultured on 60% YPD liquid medium at 30°C until the late-exponential-phase. The strain was collected by concentrated and strain density was more than 3×10^9 CFU/mL. Stressed *Z. rouxii* B-WHX-12-54 cells were washed three times with PBS by centrifugation and resuspension in PBS and then were inactivated by sterilisation. Sterilised strain and Freund's complete adjuvant were mixed (1:1) and fully emulsified as antigens. Two New Zealand White Rabbit were inoculated by subcutaneous intravenous injection with 1 mL antigens (5×10^8 CFU/mL), followed by injection of 2 mL on 5 further occasions at 10 day intervals. Then, the serum of the inoculated rabbits was collected and purified using protein A. The obtained anti-*Z. rouxii* polyclonal antibodies (pAbs) were stored at -20°C (Wang et al., 2012).

The concentration of obtained pAbs was measured using the BCA protein quantification kit (Rockford, IL, USA). Using B-WHX-12-54 as antigen, the titer of pAbs was tested by ELISA.

2.4. Western blot analysis

For the extraction of yeast proteins, the method was according to (Curto et al., 2010; Guo et al., 2016) without the acetone purification step. For bacteria proteins, 0.1 g cell pellets were added with 1 mL PBS, and then performed on an ultrasonic cell disruption system at 4°C for 30 minutes. All protein samples were quantified by BCA protein quantification kit. About 12 µL (8 µg) protein samples then mixed with 3 µL 5× SDS-PAGE sample buffers and boiled for 5 min and separated on 10% polyacrylamide gel. The separated proteins were transferred onto polyvinylidene (PVDF) membranes, the PVDF membrane was blocked for 2 h with 5% milk powder in TBST buffer. The pAbs (6.3 mg/mL) in a dilution of 1:1000 (diluted in TBST) reacted with the immobilised proteins overnight at 4°C. Plates were then washed four times with TBST, and GAR-HRP (1:2000 diluted in TBST) was added. After 2.5 h, the PVDF members were washed and developed by chemiluminescence using a ECL western blotting detection kit (Rockford, IL, USA), then recorded using a ChemiDoc XRS System (Bio-Rad Laboratories. Inc. Hercules, CA, USA) and quantified using the software Quantity One 4.6.2 software (Wang et al., 2012).

2.5. Immunofluorescence studies

The method of immunolabeling of yeast was according to Beaussart et al. (2012) and Coleman et al. (2009) with little modification. Yeast cells were washed twice by

centrifugation in PBS buffer. Cells (10^6 CFU/mL) were resuspended in 15 μ g/mL normal goat serum for 15 min at room temperature to block nonspecific antibody binding and washed twice in PBS. They were further incubated in a 18 μ g/mL solution of anti-*Z. rouxii* antibodies for 30 min at 37°C and then in a 3 μ g/mL fluorescein isothiocyanate FITC-AffiniPure Goat Anti-Rabbit IgG (H + L) for 30 min at 37°C. Images were obtained using an Andor CSU-W confocal microscopy. An 100 \times Oil objective with a DU-888U3-CSO was used for all images that were captured with Andor iQ3 live cell imaging software. The treatment of flow cytometry samples was the same as that of immunostaining samples. Flow cytometry analysis was performed with a 488-nm laser, and data were analysed using flowjo 8.5 software.

2.6. Established indirect ELISA to the detection of *Z. rouxii*

To develop a sensitive and specific ELISA, assay conditions must be optimised, including the selection of the blocking reagent (0.5% BSA, 1% BSA, 5% BSA and 5% skim milk powder in PBS), the first antibody (pAbs from Rabbit1 or Rabbit2) and the second GAR-HRP (diluted 1:1000, 1:2000, 1:3000, 1:4000 and 1:5000). The incubation times were according to Wang et al. (2012) and Li et al. (2013) reported.

The ELISA procedure used was as follows: 96-well microtiter plates (Nunc, Roskilde, Denmark) were coated with 100 μ L antigens (target or non-target strain) in PBS (pH 7.4) for 1 h at 37 °C. After incubation, plates were washed three times with 200 μ L 10 mM (PBS) with 0.05% Tween-20. Then, the plates were incubated with blocking buffer (5% skim milk in 10 mM PBS) at 37 °C for 2 h with 200 μ L per well. After three washes, the wells were blotted on absorbent paper. 100 μ L of diluted pAbs (1:4000 in PBS, 1.575 μ g/mL) were subsequently added. Following additional incubation for 1 h at 37 °C, the plate was washed three times again. Goat anti-rabbit IgG (H + L)-HRP was then added (1:2000 in PBS; 100 μ L per well) and the plate was incubated at 37°C for 1 h. After three washes, substrate solution-A TMB (1 mg/mL) and substrate solution-B (0.03% H₂O₂) were mixed (1:1) and then added to each well (150 μ L) to allow colour development for 15 min under dark conditions at 37°C. The reaction was stopped by adding 2 M of sulphuric acid (50 μ L). The OD value of each well was measured at 450 nm by a microtiter plate reader (Bio-Rad Laboratories. Inc. Hercules. CA. USA).

2.7. Cross-reaction and detection limit

Five strains of *Z. rouxii* (ATCC 2623, ATCC 8383, ATCC 66351, B-WHX-12-53 and ATCC 38531) were used for inclusivity detection. Further on, three strains of *H. uvarum* (B-WHX-12-63, B-WHX-12-14 and B-WHX-12-61), a single strain of *Z. bailli* (ATCC 8766), a single strain of *P. kluyveri* (B-WHX-12-16), a single strain of *C. orthopsilosis* (B-WHX-12-30), a single strain of *K. marxianus* (B-WHX-12-21), a single strain of *P. kudriavzevii* (B-WHX-12-36), a single strain of *C. tropicalis* (B-WHX-12-52), a single strain of *D. hansenii* (ACCC 20310), a single strain of *W. anomalus* (B-WHX-12-11), a single strain of *E. coli* (ATCC 25922) and a single strain of *A. acidoterrestris* (DSM 3923) were used as exclusivity detection strains to test the potential cross-reactivity.

ELISA signal intensity value of *Z. rouxii* B-WHX-12-54 thrice as high as background values was considered to or be regarded to be detection limit.

2.8. Detection of *Z. rouxii* in food with different sugar concentration

In order to investigate the practical applicability of the Indirect ELISA, three *Z. rouxii* (B-WHX-12-54, ATCC 2623 and ATCC 8383) artificially contaminated three food matrices (wine, apple juice concentrate and honey).

Three strains were first activated in YPD medium. 10 mL of low sugar wines were inoculated with 100 CFU/mL of the three strains. Simultaneously, three activated strains were transferred to 60% YPD medium (1% yeast extract, 2% peptone and 60% glucose) for the growth at 30 °C until the late-exponential-phase. 10 mL of high sugar of concentrated apple juice or honey were inoculated with 100 CFU/mL of the three strains. The obtained contaminated food samples were diluted 10 fold using YPD medium and cultured at 28 °C for enrichment. Then samples were taken 0, 14 and 24 h after inoculation and analysed by the developed indirect ELISA and conventional plating methods. The test was repeated three times.

3. Results and discussion

3.1. Characterisation of pAbs

The concentration of obtained pAbs was determined to be 6.1 mg/mL of R1 (Rabbit1) and 6.3 mg/mL of R2 (Rabbit2). As indicated in Supplementary Figure S1, the titer of the pAbs was 4000, resulting OD₄₅₀ values > 1. R2 is the better pAbs with high ELISA value. Therefore, pAbs of R2 was used for the establishment of an indirect ELISA system. The reactivity of the pAbs with *Z. rouxii* in different sugar stress was further confirmed by immunofluorescence studies and flow cytometry analysis (Figure 1). Both stressed-*Z. rouxii* (grown in 60% YPD) and unstressed-*Z. rouxii* (grown in YPD) stained by the pAbs exhibited distinct fluorescence (Figure 1(A)). Similar to the results in Figure 1(A), flow cytometry showed fluorescence on stressed-*Z. rouxii* and unstressed-*Z. rouxii* when the strains stained by the pAbs (Figure 1(B)). No fluorescence was observed in control sample without pAbs labelling. These results suggested that the reactivity of pAbs with *Z. rouxii*.

3.2. Establishment of indirect ELISA

Optimum results could be achieved when applying antiserum of R₂ and the HRP labeled antibodies were used in a dilution of 1:4000 and 1:2000 in PBS, respectively. The optimum condition of blocking buffer was 5% milk powder in PBS. Based on these conditions, the ELISA system was very stable (repeat several times) and the positive samples at 10⁶ CFU/mL showed ideal absorbance values higher than 1.0 and PBS blank sample is lower than 0.2.

3.3. Sensitivity, specificity, and reliability

As shown in Figure 2(A), the detection limits of the indirect ELISA for *Z. rouxii* B-WHX-12-54 was 5 × 10⁴ CFU/mL, which is comparable to the previously reported ELISA methods for the detection of *A. acidoterrestris* (10⁵ CFU/mL) (Li et al., 2013; Wang et al., 2012), *E. coli* (10⁷ CFU/mL) (Wang et al., 2007) as well as *Staphylococcus aureus* (10⁴ CFU/mL) (Pastells et al., 2015).




Figure 1. (A) Immunofluorescence images of *Z. rouxii* B-WHX-12-54 grown at different sugar concentrations. Fluorescence (b, d) and phase-contrast (a, c) images of *Z. rouxii* BW-12-54 grown in 2% and 60% YPD media, labeled with anti-*Z. rouxii* pAbs, followed FITC-conjugated secondary antibodies. (B) Histograms of negative control *Z. rouxii* B-WHX-12-54 (a), stressed-*Z. rouxii* B-WHX-12-54 (b) and unstressed-cells *Z. rouxii* B-WHX-12-54 (c) labeled with anti-*Z. rouxii* pAbs and a FITC-conjugated secondary Ab.

Few studies reported ELISA method to detect fungal, which may due to fungal share conserved antigens (Thornton & Wills, 2015). Therefore, polyclonal antibodies-based immunoassays are suitable to detect bacteria that are usually cross-reactive in the fungal detection (Schmechel et al., 2008; Vojdani, 2004). Therefore, the specificity of pAbs is particularly important to achieve ELISA detection of fungal. The specificity of the method was evaluated using 19 strains of yeast, a single strain of *E. coli*, and a single strain of *Alicyclobacillus*, part of the results are shown in Figures 2 and 3.

As shown in Figure 2, the test strain density ranged from 10^2 to 2.7×10^7 CFU/mL. When the strain density reaches the detection limit, *Z. rouxii* B-WHX-12-54 was positive. In the experimental concentration range, there was no obvious cross-reaction between the antibody and the reference strains (only *Z. bailli* ATCC 8766, *S. cerevisiae* ATCC 38531, *A. acidoterrestris* DSM 3924 and *E. coli* ATCC 25922 were shown in the Figure 2(A)). The reactivity between pAbs and the corresponding strain extracts in Figure 2(B) was also analysed by western blot. These data may explain the different affinity of pAbs with stressed-*Z. rouxii* and unstressed-*Z. rouxi*. The pAbs reacted to multiple bands in stressed-*Z. rouxii* extract and two bands in unstressed *Z. rouxii* and to two bands in *S. cerevisiae* extracts. Interestingly, *S. cerevisiae* that was weakly or nonreactive in the ELISA was cross-reactive in the western blot. The western blot uses significantly more antigen than the ELISA, and this could show reactivity to minor or underrepresented antigens. In the western blot we observed multiple antigens, some of which are to be common, but we did not specifically identify these molecules.

Moreover, the specificity of the pAbs was investigated by inclusivity detection and exclusivity evolution of the indirect ELISA system. All of three *Z. rouxii* type strains




Figure 2. (A) ELISA detection of *Z. rouxii* B-WHX-12-54 and cross-reacting analytes. Data represent the mean of three replications \pm SD. (B) Western blot analysis with the extracts of *Z. rouxii* B-WHX-12-54 and cross-reacting analytes.

(ATCC 2623 isolated from grape must and ATCC 8383 and ATCC 66351 isolated from honey) and two *Z. rouxii* isolates (B-WHX-12-53 and B-WHX-12-55 isolated from apple juice processing) not included in Figure 2 have positive results ($OD_{450} > 1$) in the test, whereas no reactivity or weakly reactivity ($OD_{450} < 0.4$) with *P. kluyveri* (B-WHX-




Figure 3. Representative images of immunostaining of fungus using an anti-*Z. rouxii* pAbs and a FITC-conjugated secondary Ab.

12-16), *K. marxianus* (B-WHX-12-21), *P. kudriavzevii* (B-WHX-12-36), *D. hansenii* (Accc 20310), *C. tropicalis* (B-WHX-12-52), *W. anomalus* (B-WHX-12-11), 2 strains of *S. cerevisiae* (1917 and 1845) and 2 strains of *H. uvarum* (B-WHX-12-14 and B-WHX-12-61) were observed. The fluorescence images of part strains are shown in [Figure 3](#), only *C. orthopsilosis* (B-WHX-12-30) showed weakly fluorescence. The immunofluorescence studies confirmed the ELISA results. Most of the tested strains used in this study are isolated from apple orchards and apple juice processing which have the potential for spoilage (Wang et al., 2015a). In total, these results suggested that the pAbs-based ELISA can recognise the common food fungus contamination, and *Z. rouxii* contribute the most important to ELISA signal. Thus, pAs-based ELISA for detection of *Z. rouxii* has high sensitivity and specificity and has application value.

3.4. Detection of *Z. rouxii* in food with different sugar conditions

A preliminary measurement of the ELISA for the direct detection of *Z. rouxii* in concentrated apple juice revealed that high density of *Z. rouxii* in high soluble solids content (SSC) samples resulted in lower absorbance values. Therefore, the SSC of the sample is an important factor affecting the ELISA results. The result may be due to the *Z. rouxii* cells in the high SSC concentrated apple juice can't bind to the bottom of microtiter plates due to the strong viscosity of the juice. The spoilage of *Z. rouxii* usually found on the surface of concentrated apple juice (Wang et al., 2015a) also supports this result. This issue will be solved since the SSC of diluted food samples changed ranging from (4 °Brix to 75 °Brix) to (5 °Brix to 15 °Brix) (Supplementary Table S1), which is suitable for ELISA detection.

An indirect ELISA developed in food (concentrated apple juice, wine, and honey) was evaluated for the detection of *Z. rouxii* in contaminated samples. In preliminary tests, the food samples (uncontaminated or spiked with 10^7 CFU/mL *Z. rouxii*) were serially diluted in PBS and analysed by the indirect ELISA. Detection limits for *Z. rouxii* were in the range of 10^5 CFU/mL, thus confirming the ability of the indirect ELISA method for the sensitive detection of *Z. rouxii* in food samples. However, the background absorbance values measured for the uncontaminated food samples varied depending on the kinds of the foods ([Table 2](#)). Analysis of wine and concentrated apple juice resulted in low background values of about 0.07-0.19, while the honey exhibited higher background values of about 0.53-0.63. One reason could be that the composition in honey should be consideration to interfere the results. However, the background values of honey could be diminished to proper values by a simple dilution 4 fold in PBS. Therefore, the dilution step results in a decrease of target concentration, which could be solved by an extended *Z. rouxii* enrichment time of honey samples.

In practical ELISA, for each tested strains, samples after contamination for 0, 14 and 24 h were simultaneously analysed by ELISA system and conventional plating methods. As expected, concentrated apple juice and wine samples could be detected within 24 h ([Table 2](#)). At this time, the plating results showed that the strain concentration was ranging from about 5×10^4 CFU/mL to 2×10^5 CFU/mL (Supplementary Table S2) and the ELISA signal intensity of food samples was at least higher than 3 fold of the samples at 0 h. While for Honey, samples could be detected positive within 30 h as indicated in [Table 3](#). The developed indirect ELISA possesses the capacity of detecting

Table 2. ELISA results for different drinks (diluted 10 fold in YPD medium) contaminated artificially with about 100 CFU/mL of isolate B-WHX-12-54 or ATCC 2623 and ATCC 8383.

	Apple juice concentrate			Wine			Honey		
	B-WHX-12-54	ATCC2623	ATCC8383	B-WHX-12-53	ATCC2623	ATCC8383	B-WHX-12-53	ATCC2623	ATCC8383
0 h	0.1 ± 0.01	0.07 ± 0.01	0.08 ± 0.01	0.19 ± 0.01	0.14 ± 0.03	0.12 ± 0.01	0.68 ± 0.02	0.53 ± 0.07	0.63 ± 0.06
14 h	0.11 ± 0.01	0.09 ± 0.01	0.13 ± 0.02	0.21 ± 0.05	0.19 ± 0.02	0.16 ± 0.02	0.85 ± 0.07	0.86 ± 0.08	0.71 ± 0.07
24 h	1.09 ± 0.03	0.4 ± 0.04	0.46 ± 0.02	0.51 ± 0.05	0.4 ± 0.02	0.36 ± 0.02	1.56 ± 0.07	1.77 ± 0.08	1.37 ± 0.08

Table 3. ELISA results for honey (diluted 10 fold in YPD medium) contaminated artificially with about 100 CFU/mL of isolate B-WHX-12-54 or ATCC 2623 and ATCC 8383.

Dilution 4 fold in PBS	Honey		
	B-WHX-12-53	ATCC2623	ATCC8383
0 h	0.151 ± 0.04	0.1 ± 0.02	0.11 ± 0.02
14 h	0.14 ± 0.05	0.09 ± 0.01	0.11 ± 0.02
24 h	0.18 ± 0.02	0.10 ± 0.01	0.12 ± 0.02
30 h	1.72 ± 0.09	1.76 ± 0.01	1.56 ± 0.01

Z. rouxii directly from contaminated food samples after a 24 h or 30 h incubation step. As low as 100 CFU/mL of *Z. rouxii* can be detected. The strain density is equal to previously reported e-nose combined with chemometrics method for the detection of *Z. rouxii* (10^2 CFU/mL). But the e-nose method is limited to predicting the concentration of *Z. rouxii* in concentrated apple juice. Besides, the method requires that the sample pretreatment is not suitable for routine testing. It is worth noting that detection of *Z. rouxii* in high sugar selective medium requires a very long time (96-120 h). The established indirect ELISA here only requires 24-30 h enrichment time and simple to operate. These results demonstrated established indirect ELISA can be used to rapid detect *Z. rouxii* in food samples.

4. Conclusions

We observed the anti-*Z. rouxii* pAbs could be obviously recognise *Z. rouxii* and did not or weakly react to spoilage yeast or bacteria in different food. Here in, we reported a highly sensitive and reliable polyclonal antibodies-based indirect ELISA for *Z. rouxii* detection in different food. It could be demonstrating that pAbs recognise *Z. rouxii* growth in different food containing different sugar stress condition, probably due to the conserved expression of *Z. rouxii* surface proteins. As far as we know, the detection of osmotolerant yeast depending on the high sugar selective medium is a long time-consuming classical microbiological method. Due to the high affinity of Abs and the ability to bind with various nanomaterials, the antibodies-based bioanalytical techniques could be used to ultimately replace classical methods to ensure food quality.

Acknowledgements

The authors would like to acknowledge Mrs Wang, who is responsible for the Life science research platform laser confocal chamber of Northwest A & F University. She generously provided the microscopy technology.

Disclosure statement

No potential conflict of interest was reported by the authors.

Funding

This work was supported by National Natural Science Foundation of China [31371814]; major projects of strategic emerging industries in Shaanxi province of China [2016KTCQ03-12].

References

- Arroyomanzanas, N., Huertaspérez, J. F., Gámizgracia, L., & Garcíacampaña, A. M. (2015). Simple and efficient methodology to determine mycotoxins in cereal syrups. *Food Chemistry*, *177*, 274–279.
- Beaussart, A., Alsteens, D., El-Kirat-Chatel, S., Lipke, P. N., Kucharíková, S., Van Dijck, P., & Dufrière, Y. F. (2012). Single-molecule imaging and functional analysis of Als adhesins and mannans during *Candida albicans* morphogenesis. *ACS Nano*, *6*, 10950–10964.
- Coleman, D. A., Oh, S. H., Zhao, X. M., Zhao, H. Y., Hutchins, J. T., Vernachio, J. H., ... Hoyer, L. L. (2009). Monoclonal antibodies specific for *Candida albicans* Als3 that immunolabel fungal cells in vitro and in vivo and block adhesion to host surfaces. *Journal of Microbiological Methods*, *78*, 71–78.
- Curto, M., Valledor, L., Navarrete, C., Gutierrez, D., Sychrova, H., Ramos, J., & Jorrián, J. (2010). 2-DE based proteomic analysis of *Saccharomyces cerevisiae* wild and K⁺ transport-affected mutant (trk1,2) strains at the growth exponential and stationary phases. *Journal of Proteomics*, *73*, 2316–2335.
- Dakal, T. C., Solieri, L., & Giudici, P. (2014). Adaptive response and tolerance to sugar and salt stress in the food yeast *Zygosaccharomyces rouxii*. *International Journal of Food Microbiology*, *185*, 140–157.
- Gobbi, E., Falasconi, M., Concina, I., Mantero, G., Bianchi, F., Mattarozzi, M., & Sberveglieri, G. (2010). Electronic nose and *Alicyclobacillus* spp. spoilage of fruit juices: an emerging diagnostic tool. *Food Control*, *21*, 374–382.
- Gobbi, E., Falasconi, M., Zambotti, G., Sberveglieri, V., Pulvirenti, A., & Sberveglieri, G. (2015). Rapid diagnosis of Enterobacteriaceae in vegetable soups by a metal oxide sensor based electronic nose. *Sensors and Actuators B: Chemical*, *207*, 1104–1113.
- Guo, H., Niu, C., Liu, B., Wei, J. P., Wang, H. X., Yuan, Y. H., & Yue, T. L. (2016). Protein abundance changes of *Zygosaccharomyces rouxii* in different sugar concentrations. *International Journal of Food Microbiology*, *233*, 44–51.
- Guo, H., Yuan, Y. H., Niu, C., Wang, Z. L., Qiu, Y., & Yue, T. L. (2017). Wash-free colorimetric homogeneous immunoassay for *Zygosaccharomyces rouxii*. *RSC Advances*, *7*, 34307–34314.
- Harrison, E., Muir, A., Stratford, M., & Wheals, A. (2011). Species-specific PCR primers for the rapid identification of yeasts of the genus *Zygosaccharomyces*. *FEMS Yeast Research*, *11*, 356–365.
- Hoffmann, H., Baldofski, S., Hoffmann, K., Flemig, S., Silva, C. P., Esteves, V. I., & Schneider, R. J. (2016). Structural considerations on the selectivity of an immunoassay for sulfamethoxazole. *Talanta*, *158*, 198–207.
- Hulin, M., & Wheals, A. (2014). Rapid identification of *Zygosaccharomyces* with genus-specific primers. *International Journal of Food Microbiology*, *173*, 9–13.
- Kumar, B. K., Raghunath, P., Devegowda, D., Deekshit, V. K., Venugopal, M. N., & Karunasagar, I. (2011). Development of monoclonal antibody based sandwich ELISA for the rapid detection of pathogenic *Vibrio parahaemolyticus* in seafood. *International Journal of Food Microbiology*, *145*, 244–249.
- Li, J. K., Huang, R., Xia, K., & Liu, L. (2014). Double antibodies sandwich enzyme-linked immunosorbent assay for the detection of *Alicyclobacillus acidoterrestris* in apple juice concentrate. *Food Control*, *40*, 172–176.
- Li, J. K., Xia, K., & Yu, C. Z. (2013). Detection of *Alicyclobacillus acidoterrestris* in apple juice concentrate by enzyme-linked immunosorbent assay. *Food Control*, *30*, 251–254.
- Loureiro, V., & Malfeitoferreira, M. (2003). Spoilage yeasts in the wine industry. *International Journal of Food Microbiology*, *86*, 23–50.
- Martorell, P., Stratford, M., Steels, H., Fernandez-Espinar, M. T., & Querol, A. (2007). Physiological characterization of spoilage strains of *Zygosaccharomyces bailii* and *Zygosaccharomyces rouxii* isolated from high sugar environments. *International Journal of Food Microbiology*, *114*, 234–242.
- Mast, S., Dietrich, R., Didier, A., & Martlbauer, E. (2016). Development of a polyclonal antibody-based sandwich enzyme-linked immunosorbent assay for the detection of spores of

- alicyclobacillus acidoterrestris in various fruit juices. *Journal of Agriculture and Food Chemistry*, 64, 497–504.
- Pastells, C., Acosta, G., Pascual, N., Albericio, F., Royo, M., & Marco, M. P. (2015). An immunochemical strategy based on peptidoglycan synthetic peptide epitopes to diagnose *Staphylococcus aureus* infections. *Analytica Chimica Acta*, 889, 203–211.
- Pearson, B., & McKee, R. (1992). Rapid identification of *Saccharomyces cerevisiae*, *Zygosaccharomyces bailii* and *Zygosaccharomyces rouxii*. *International Journal of Food Microbiology*, 16, 63–67.
- Schmechel, D., Green, B. J., Blachere, F. M., Janotka, E., & Beezhold, D. H. (2008). Analytical bias of cross-reactive polyclonal antibodies for environmental immunoassays of *Alternaria alternata*. *Journal of Allergy Clinical Immunology*, 121, 763–768.
- Thornton, C. R., & Wills, O. E. (2015). Immunodetection of fungal and oomycete pathogens: Established and emerging threats to human health, animal welfare and global food security. *Critical Reviews in Microbiology*, 41, 27–51.
- Vojdani, A. (2004). Cross-reactivity of *Aspergillus*, *Penicillium*, and *Stachybotrys* antigens using affinity-purified antibodies and immunoassay. *International Archives Occupational Environmental Health*, 59, 256–265.
- Wang, H. X., Hu, Z. Q., Long, F. Y., Guo, C. F., Yuan, Y. H., & Yue, T. L. (2015). Detection of *Zygosaccharomyces rouxii* and *Candida tropicalis* in a high-sugar medium by a metal oxide sensor-based electronic nose and comparison with test panel evaluation. *Journal of Food Protection*, 78, 2052–2063.
- Wang, H. X., Hu, Z. Q., Long, F. Y., Niu, C., Guo, C. F., Yuan, Y. H., & Yue, T. L. (2015). Characterization of osmotolerant yeasts and yeast-like molds from apple orchards and apple juice processing plants in China and investigation of their spoilage potential. *Journal of Food Science*, 80, M1850–M1860.
- Wang, H. X., Hu, Z. Q., Long, F. Y., Guo, C. F., Yuan, Y. H., & Yue, T. L. (2016). Early detection of *Zygosaccharomyces rouxii* – spawned spoilage in apple juice by electronic nose combined with chemometrics. *International Journal of Food Microbiology*, 217, 68–78.
- Wang, N., He, M., & Shi, H. C. (2007). Novel indirect enzyme-linked immunosorbent assay (ELISA) method to detect total *E. coli* in water environment. *Analytica Chimica Acta*, 590, 224–231.
- Wang, Z., Beier, R. C., Sheng, Y., Zhang, S., Jiang, W., Wang, Z., & Shen, J. (2013). Monoclonal antibodies with group specificity toward sulfonamides: Selection of hapten and antibody selectivity. *Analytical Bioanalytical Chemistry*, 405, 4027–4037.
- Wang, Z. L., Yue, T. L., Yuan, Y. H., Cai, R., Guo, C. F., Wang, X., & Niu, C. (2012). Development of polyclonal antibody-based indirect enzyme-linked immunosorbent assay for the detection of *Alicyclobacillus* strains in apple juice. *Journal of Food Science*, 77, M643–M649.
- Watanabe, J., Uehara, K., & Mogi, Y. (2013). Adaptation of the osmotolerant yeast *Zygosaccharomyces rouxii* to an osmotic environment through copy number amplification of FLO11D. *Genetics*, 195, 393–405.
- Zhang, J. B., Yue, T. L., & Yuan, Y. H. (2013). *Alicyclobacillus* contamination in the production line of kiwi products in China. *PloS One*, 8, e67704.
- Zuehlke, J. M., Petrova, B., & Edwards, C. G. (2013). Advances in the control of wine spoilage by *Zygosaccharomyces* and *Dekkera/Brettanomyces*. *Annual Review of Food Science and Technology*, 4, 57–78.